

Estd. 1966

®



DELIVERING TO THE WORLD

# 12 Channel Fully Autoclavable Micropipette (MODEL NO : MPMC 301)



**12 Channel Fully Autoclavable Micropipette Manufactured by Desco Medical India**



Decide with Confidence

## 12 Channel Fully Autoclavable Micropipette

Estd. 1966

®



DELIVERING TO THE WORLD



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

### Specifications:-

This is Most accurate, reliable, comfortable for extensive pipetting in today's research. Research today demands only highest quality when it comes to pipetting and precision. Handypipette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.



# 12 Channel Fully Autoclavable Micropipette (MODEL NO : MPMC 301)



DELIVERING TO THE WORLD

- \* Lower assembly rotates 360 for optimal working position.
- \* Optimum piston-spring balance for soft plunger movement.
- \* Easy Volume setting.
- \* Easy sequential tip ejection system.
- \* Easy user calibration.

We are Manufacturer Exporters and Suppliers of 12 Channel Fully Autoclavable Micropipette at wholesale competitive prices from India

<b>MODEL NO:</b>	<b>VOLUME RANGE</b>	<b>ACCURACY %</b>	<b>PRECISION(CV%)</b>	<b>PACK SIZE</b>
<b>MPMC 301</b>	0.5 - 10 ul	4.00	4.00	1 Pc.
		2.50	2.50	
		1.50	1.50	
<b>MPMC 302</b>	2.0 - 20 ul	4.00	4.00	1 Pc.
		1.50	1.00	
		1.50	1.00	
<b>MPMC 303</b>	5.0 - 50 ul	3.00	2.00	1 Pc.
		1.50	1.00	
		0.50	0.50	
<b>MPMC 304</b>	10 - 100 ul	3.00	1.50	1 Pc.
		1.00	0.50	
		1.00	0.50	
<b>MPMC 305</b>	20 - 200 ul	2.00	0.80	1 Pc.
		1.00	0.50	
		0.70	0.25	
<b>MPMC 306</b>	30 - 300 ul	2.00	0.80	1 Pc.
		1.00	0.50	
		0.70	0.25	