

Estd. 1966

®



DELIVERING TO THE WORLD

8 Channel Fully Autoclavable Micropipette (MODEL NO : MPMC 101)



8 Channel Fully Autoclavable Micropipette Manufactured by Desco Medical India



Decide with Confidence

8 Channel Fully Autoclavable Micropipette

Estd. 1966

®



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

This is Most accurate, reliable, comfortable for extensive pipetting in today's research. Research today demands only highest quality when it comes to pipetting and precision. Handypipette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

Estd. 1966



DELIVERING TO THE WORLD

8 Channel Fully Autoclavable Micropipette (MODEL NO : MPMC 101)



- * Lower assembly rotates 360 for optimal working position.
- * Optimum piston-spring balance for soft plunger movement.
- * Easy Volume setting.
- * Easy sequential tip ejection system.
- * Easy user calibration.

We are Manufacturer Exporters and Suppliers of 8 Channel Fully Autoclavable Micropipette at wholesale competitive prices from India

MODEL NO:	VOLUME RANGE	ACCURACY %	PRECISION(CV%)	PACK SIZE
MPMC 101	0.5 - 10 ul	4.00	4.00	1 Pc.
		2.50	2.50	
		1.50	1.50	
MPMC 102	2.0 - 20 ul	4.00	4.00	1 Pc.
		1.50	1.00	
		1.50	1.00	
MPMC 103	5.0 - 50 ul	3.00	2.00	1 Pc.
		1.50	1.00	
		0.50	0.50	
MPMC 104	10 - 100 ul	3.00	1.50	1 Pc.
		1.00	0.50	
		1.00	0.50	
MPMC 105	20 - 200 ul	2.00	0.80	1 Pc.
		1.00	0.50	
		0.70	0.25	
MPMC 106	30 - 300 ul	2.00	0.80	1 Pc.
		1.00	0.50	
		0.70	0.25	