

Estd. 1966

®



DELIVERING TO THE WORLD

8 Channel Semi Autoclavable Micropipette (MODEL NO : MPMC 201)



8 Channel Semi Autoclavable Micropipette Manufactured by Desco Medical India



Decide with Confidence

8 Channel Semi Autoclavable Micropipette

Estd. 1966

®



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

This is Most accurate, reliable, comfortable for extensive pipetting in today's research. Research today demands only highest quality when it comes to pipetting and precision. Handypipette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)



8 Channel Semi Autoclavable Micropipette (MODEL NO : MPMC 201)



- * Lower assembly rotates 360 for optimal working position.
- * Optimum piston-spring balance for soft plunger movement.
- * Easy Volume setting.
- * Easy sequential tip ejection system.
- * Easy user calibration.

We are Manufacturer Exporters and Suppliers of 8 Channel Semi Autoclavable Micropipette at wholesale competitive prices from India

MODEL NO:	VOLUME RANGE	ACCURACY %	PRECISION(CV%)	PACK SIZE
MPMC 201	0.5 - 10 ul	4.00	4.00	1 Pc.
		2.50	2.50	
		1.50	1.50	
MPMC 202	2.0 - 20 ul	4.00	4.00	1 Pc.
		1.50	1.00	
		1.50	1.00	
MPMC 203	5.0 - 50 ul	3.00	2.00	1 Pc.
		1.50	1.00	
		0.50	0.50	
MPMC 204	10 - 100 ul	3.00	1.50	1 Pc.
		1.00	0.50	
		1.00	0.50	
MPMC 205	20 - 200 ul	2.00	0.80	1 Pc.
		1.00	0.50	
		0.70	0.25	
MPMC 206	30 - 300 ul	2.00	0.80	1 Pc.
		1.00	0.50	
		0.70	0.25	