

Estd. 1966



DELIVERING TO THE WORLD

# Autoclave Pressure Cooker Aluminium Electric (MODEL NO : AUPC 101)



## Autoclave Pressure Cooker Electric Made of Aluminium Manufactured by Desco Medical India



Decide with Confidence

### Autoclave Pressure Cooker Aluminium Electric

Estd. 1966



DELIVERING TO THE WORLD



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.autoclavesterilizersindia.com](http://www.autoclavesterilizersindia.com)

[www.descoonline.in](http://www.descoonline.in)



# Autoclave Pressure Cooker Aluminium Electric (MODEL NO : AUPC 101)



## Autoclave Pressure Cooker Electric Made of Aluminium Manufactured by Desco Medical India

### Specifications :-

Autoclave Pressure Cooker is also known as Pressure Steam Sterilizers, Steam Autoclave, Portable Autoclaves Sterilizers

Autoclave Pressure Cooker is a part of Steam Sterilizing Equipment

Autoclave Pressure Cooker is Electrical and has mirror finish

The body and lid has Seamless Construction from Extra Thick Aluminium Sheet Body & Lid for more safety and life long performance.

The Anti-Bulging & Fuel Saver Base is provided for Additional Strength, Fuel Saving and Time Saving.

The depression at centre of the lid helps in Additional Strength.

The Gasket is ISI Marked food grade synthetic rubber gasket for better reliability, hygiene & longer life.

The Gasket Release System (GRS) & fusible plug type safety valve is available for Triple-Safety.

The pressure gauge is Colour-coded Pressure Gauge to indicate the safe and danger zones.

The repeating type safety valve is provided for better Safety and longer & error-free performance.

The vacuum release & steam release valves are made of Stainless Steel for longer life and to avoid toxic deposits on valves due to oxidation.

The heating element is ISI Marked which is made of thick and heavy copper pipe with very thick chrome plating.

The Electric Cable has heavy electric wires to comfortably bear the heating load.

The 3 Pin Plug & Element Connector is of High quality for electrical safety.

The handles is made from heat resistant, impact tolerant, finest quality bakelite, fitted with Stainless Steel screws to avoid easy loosening of handles.

The inside tripod stand is made of Stainless Steel with extra height to avoid wetting of instruments.

The outer stand is well designed to provide proper stability to the Autoclave Pressure Cooker.

We are Manufacturer, Exporters and Suppliers of Autoclave Pressure Cooker Aluminium Electric at wholesale competitive prices from India.

Model No	Size in inch	Size in mm	Capacity (Ltr)	Heating Element
AUPC 101	9 x 9 Inch	225 mm x 225 mm	10 Ltr	1 KW
AUPC 102	9 x 11 Inch	225 mm x 275 mm	13 Ltr	1 KW
AUPC 103	9 x 13 Inch	225 mm x 325 mm	15 Ltr	1.5 KW
AUPC 103.1	12 x 11 Inch	300 mm x 185 mm	18 Ltr	1.5 KW
AUPC 103.2	12 x 11 Inch	300 mm x 220 mm	20 Ltr	1.5 KW
AUPC 104	12 x 12 Inch	300 mm x 300 mm	23.5 Ltr	2 KW
AUPC 105	12 x 14 Inch	300 mm x 350 mm	26 Ltr	2 KW
AUPC 106	12 x 16 Inch	300 mm x 400 mm	30 Ltr	2 KW
AUPC 107	12 x 18 Inch	300 mm x 450 mm	33.5 Ltr	2 KW
AUPC 108	14 x 15 Inch	350 mm x 375 mm	36 Ltr	2 KW
AUPC 109	12 x 22 Inch	300 mm x 550 mm	40 Ltr	2 KW
AUPC 110	14 x 17 Inch	350 mm x 425 mm	41 Ltr	2 KW
AUPC 111	14 x 19 Inch	350 mm x 475 mm	45 Ltr	2 KW
AUPC 112	14 x 21 Inch	350 mm x 525 mm	50 Ltr	2 KW
AUPC 113	14 x 24 Inch	350 mm x 600 mm	60 Ltr	2 KW