



Autoclave Pressure Cooker Aluminium Non Electric (MODEL NO : AUPC 151)



**Autoclave Pressure Cooker Non Electric Made of
Aluminium Manufactured by Desco Medical India**



Decide with Confidence

Autoclave Pressure Cooker Aluminium Non Electric



www.descomedicalindia.com

www.autoclavesterilizersindia.com

www.descoonline.in

Autoclave Pressure Cooker Non Electric Made of Aluminium Manufactured by Desco Medical India

Specifications :-

Autoclave Pressure Cooker is also known as Pressure Steam Sterilizers, Steam Autoclave, Portable Autoclaves Sterilizers
 Autoclave Pressure Cooker is a part of Steam Sterilizing Equipment
 Autoclave Pressure Cooker is Non Electric and has mirror finish
 The body and lid has Seamless Construction from Extra Thick Aluminium Sheet Body & Lid for more safety and life long performance.
 The Anti-Bulging & Fuel Saver Base is provided for Additional Strength, Fuel Saving and Time Saving.
 The depression at centre of the lid helps in Additional Strength.
 The Gasket is ISI Marked food grade synthetic rubber gasket for better reliability, hygiene & longer life.
 The Gasket Release System (GRS) & fusible plug type safety valve is available for Triple-Safety.
 The pressure gauge is Colour-coded Pressure Gauge to indicate the safe and danger zones.
 The repeating type safety valve is provided for better Safety and longer & error-free performance.
 The vacuum release & steam release valves are made of Stainless Steel for longer life and to avoid toxic deposits on valves due to oxidation.
 The handles is made from heat resistant, impact tolerant, finest quality bakelite, fitted with Stainless Steel screws to avoid easy loosening of handles.
 The inside tripod stand is made of Stainless Steel with extra height to avoid wetting of instruments.
 The outer stand is well designed to provide proper stability to the Autoclave Pressure Cooker.
 We are Manufacturer, Exporters and Suppliers of Autoclave Pressure Cooker Aluminium Non Electric at wholesale competitive prices from India.

Model No	Size in inch	Size in mm	Capacity (Ltr)
AUPC 151	9 x 9 Inch	225 mm x 225 mm	10 Ltr
AUPC 152	9 x 11 Inch	225 mm x 275 mm	13 Ltr
AUPC 153	9 x 13 Inch	225 mm x 325 mm	15 Ltr
AUPC 153.1	12 x 11 Inch	300 mm x 185 mm	18 Ltr
AUPC 153.2	12 x 11 Inch	300 mm x 220 mm	20 Ltr
AUPC 154	12 x 12 Inch	300 mm x 300 mm	23.5 Ltr
AUPC 155	12 x 14 Inch	300 mm x 350 mm	26 Ltr
AUPC 156	12 x 16 Inch	300 mm x 400 mm	30 Ltr
AUPC 157	12 x 18 Inch	300 mm x 450 mm	33.5 Ltr
AUPC 158	14 x 15 Inch	350 mm x 375 mm	36 Ltr
AUPC 159	12 x 22 Inch	300 mm x 550 mm	40 Ltr
AUPC 160	14 x 17 Inch	350 mm x 425 mm	41 Ltr
AUPC 161	14 x 19 Inch	350 mm x 475 mm	45 Ltr
AUPC 162	14 x 21 Inch	350 mm x 525 mm	50 Ltr
AUPC 163	14 x 24 Inch	350 mm x 600 mm	60 Ltr