



Autoclave Pressure Cooker Stainless Steel Electric (MODEL NO : AUPC 301)



Autoclave Pressure Cooker Electric Made of Stainless Steel Manufactured by Desco Medical India



Decide with Confidence

Autoclave Pressure Cooker Stainless Steel Electric

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.autoclavesterilizersindia.com

www.descoonline.in

Autoclave Pressure Cooker Electric Made of Stainless Steel Manufactured by Desco Medical India

Specifications :-

Autoclave Pressure Cooker P Type is also known as Pressure Steam Sterilizers, Steam Autoclave, Portable Autoclaves Sterilizers

Autoclave Pressure Cooker is a part of Steam Sterilizing Equipment

Autoclave Pressure Cooker is Electrical and is made of Stainless Steel

The body and lid has Seamless Construction from Extra Thick Stainless Steel Sheet Body & Lid for more safety and life long performance.

The Anti-Bulging & Fuel Saver Base is provided for Additional Strength, Fuel Saving and Time Saving.

The depression at centre of the lid helps in Additional Strength.

The Gasket is ISI Marked food grade synthetic rubber gasket for better reliability, hygiene & longer life.

The Gasket Release System (GRS) & fusible plug type safety valve is available for Triple-Safety.

The pressure gauge is Colour-coded Pressure Gauge to indicate the safe and danger zones.

The repeating type safety valve is provided for better Safety and longer & error-free performance.

The vacuum release & steam release valves are made of Stainless Steel for longer life and to avoid toxic deposits on valves due to oxidation.

The heating element is ISI Marked which is made of thick and heavy copper pipe with very thick chrome plating.

The Electric Cable has heavy electric wires to comfortably bear the heating load.

The 3 Pin Plug & Element Connector is of High quality for electrical safety.

The handles is made from heat resistant, impact tolerant, finest quality bakelite, fitted with Stainless Steel screws to avoid easy loosening of handles.

The inside tripod stand is made of Stainless Steel with extra height to avoid wetting of instruments.

The outer stand is well designed to provide proper stability to the Autoclave Pressure Cooker.

We are Manufacturer, Exporters and Suppliers of Autoclave Pressure Cooker Electric Made of Stainless Steel at wholesale competitive prices from India.

Model No	Size in inch	Size in mm	Capacity (Ltr)	Heating Element
AUPC 301	9 X 9 Inch	225 X 225 mm	9 Ltr	1 KW
AUPC 302	9 X 11 Inch	225 X 275 mm	11 Ltr	1 KW
AUPC 303	9 X 16 Inch	225 X 400 mm	16 Ltr	1.5 KW
AUPC 304	12 X 12 Inch	300 X 300 mm	21 Ltr	2 KW
AUPC 305	12 X 14 Inch	300 X 350 mm	25 Ltr	2 KW
AUPC 306	14 X 13 Inch	350 X 525 mm	34 Ltr	2 KW
AUPC 307	14 X 15 Inch	350 X 375 mm	39 Ltr	2 KW
AUPC 308	12 X 22 Inch	350 X 550 mm	40 Ltr	2 KW