



Autoclave Pressure Cooker Stainless Steel Non Electric (MODEL NO : AUPC 351)



Autoclave Pressure Cooker Non Electric Made of Stainless Steel Manufactured by Desco Medical India



Decide with Confidence

**Autoclave Pressure Cooker Stainless Steel
Non Electric**

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.autoclavesterilizersindia.com

www.descoonline.in



Autoclave Pressure Cooker Stainless Steel Non Electric (MODEL NO : AUPC 351)



Autoclave Pressure Cooker Non Electric Made of Stainless Steel Manufactured by Desco Medical India

Specifications :-

Autoclave Pressure Cooker P Type is also known as Pressure Steam Sterilizers, Steam Autoclave, Portable Autoclaves Sterilizers

Autoclave Pressure Cooker is a part of Steam Sterilizing Equipment

Autoclave Pressure Cooker is Non Electric and is made of Stainless Steel

The body and lid has Seamless Construction from Extra Thick Stainless Steel Sheet Body & Lid for more safety and life long performance.

The Anti-Bulging & Fuel Saver Base is provided for Additional Strength, Fuel Saving and Time Saving.

The depression at centre of the lid helps in Additional Strength.

The Gasket is ISI Marked food grade synthetic rubber gasket for better reliability, hygiene & longer life.

The Gasket Release System (GRS) & fusible plug type safety valve is available for Triple-Safety.

The pressure gauge is Colour-coded Pressure Gauge to indicate the safe and danger zones.

The repeating type safety valve is provided for better Safety and longer & error-free performance.

The vacuum release & steam release valves are made of Stainless Steel for longer life and to avoid toxic deposits on valves due to oxidation.

The handles is made from heat resistant, impact tolerant, finest quality bakelite, fitted with Stainless Steel screws to avoid easy loosening of handles.

The inside tripod stand is made of Stainless Steel with extra height to avoid wetting of instruments.

The outer stand is well designed to provide proper stability to the Autoclave Pressure Cooker.

We are Manufacturer, Exporters and Suppliers of Autoclave Pressure Cooker Non Electric Made of Stainless Steel at wholesale competitive prices from India.

Model No	Size in inch	Size in mm	Capacity (Ltr)
AUPC 351	9 X 9 Inch	225 X 225 mm	9 Ltr
AUPC 352	9 X 11 Inch	225 X 275 mm	11 Ltr
AUPC 353	9 X 16 Inch	225 X 400 mm	16 Ltr
AUPC 354	12 X 12 Inch	300 X 300 mm	21 Ltr
AUPC 355	12 X 14 Inch	300 X 350 mm	25 Ltr
AUPC 356	14 X 13 Inch	350 X 525 mm	34 Ltr
AUPC 357	14 X 15 Inch	350 X 375 mm	39 Ltr
AUPC 358	12 X 22 Inch	350 X 550 mm	40 Ltr