



Autoclave Wing Nut Aluminium Electric (MODEL NO : AUWN 101)



**Autoclave Wing Nut Electric Made of Aluminium
Manufactured by Desco Medical India**



Decide with Confidence

Autoclave Wing Nut Aluminium Electric

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.autoclavesterilizersindia.com

www.descoonline.in

Autoclave Wing Nut Electric Made of Aluminium Manufactured by Desco Medical India

Specifications :-

- Autoclave Wing Nut is also known as Wing Nut Sterilizers, Portable Wing Nut Autoclaves Sterilizers
 - Autoclave Wing Nut is a part of Steam Sterilizing Equipment
 - Autoclave Wing Nut is Electrical and has mirror finish
 - The body and lid has Seamless Construction from Extra Thick Aluminium Sheet Body & Lid for more safety and life long performance.
 - The Anti-Bulging & Fuel Saver Base is provided for Additional Strength, Fuel Saving and Time Saving.
 - The depression at centre of the lid helps in Additional Strength.
 - The Gasket is ISI Marked food grade synthetic rubber gasket for better reliability, hygiene & longer life.
 - The Gasket Release System (GRS) & fusible plug type safety valve is available for Triple-Safety.
 - The pressure gauge is Colour-coded Pressure Gauge to indicate the safe and danger zones.
 - The repeating type safety valve is provided for better Safety and longer & error-free performance.
 - The vacuum release & steam release valves are made of Stainless Steel for longer life and to avoid toxic deposits on valves due to oxidation.
 - The heating element is ISI Marked which is made of thick and heavy copper pipe with very thick chrome plating.
 - The Electric Cable has heavy electric wires to comfortably bear the heating load.
 - The 3 Pin Plug & Element Connector is of High quality for electrical safety.
 - The handles is made from heat resistant, impact tolerant, finest quality bakelite, fitted with Stainless Steel screws to avoid easy loosening of handles.
 - The inside tripod stand is made of Stainless Steel with extra height to avoid wetting of instruments.
 - The outer stand is well designed to provide proper stability to the Autoclave Wing Nut.
- We are Manufacturer, Exporters and Suppliers of Autoclave Wing Nut Electric Made of Aluminium at wholesale competitive prices from India.

Model No	Size in inch	Size in mm	Capacity (Ltr)	Heating Element
AUWN 101	9 X 9 Inch	230 X 230 mm	10 Ltr	2 KW
AUWN 102	9 X 13 Inch	230 X 330 mm	15 Ltr	2 KW
AUWN 103	12 X 12 Inch	300 X 300 mm	20 Ltr	2 KW
AUWN 104	12 X 12 Inch	300 X 300 mm	21 Ltr	2 KW
AUWN 105	12 X 14 Inch	300 X 350 mm	25 Ltr	2 KW
AUWN 106	12 X 15 Inch	300 X 375 mm	27 Ltr	2 KW
AUWN 107	12 X 16 Inch	300 X 400 mm	30 Ltr	2 KW
AUWN 108	14 X 15 Inch	350 X 375 mm	39 Ltr	2 KW