



# Autoclave Wing Nut Aluminium Non Electric (MODEL NO : AUWN 151)



**Autoclave Wing Nut Non Electric Made of Aluminium  
Manufactured by Desco Medical India**



Decide with Confidence

**Autoclave Wing Nut Aluminium Non Electric**



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.autoclavesterilizersindia.com](http://www.autoclavesterilizersindia.com)

[www.descoonline.in](http://www.descoonline.in)

## Autoclave Wing Nut Non Electric Made of Aluminium Manufactured by Desco Medical India

### Specifications :-

- Autoclave Wing Nut is also known as Wing Nut Sterilizers, Portable Wing Nut Autoclaves Sterilizers
- Autoclave Wing Nut is a part of Steam Sterilizing Equipment
- Autoclave Wing Nut is Non Electric and has mirror finish
- The body and lid has Seamless Construction from Extra Thick Aluminium Sheet Body & Lid for more safety and life long performance.
- The Anti-Bulging & Fuel Saver Base is provided for Additional Strength, Fuel Saving and Time Saving.
- The depression at centre of the lid helps in Additional Strength.
- The Gasket is ISI Marked food grade synthetic rubber gasket for better reliability, hygiene & longer life.
- The Gasket Release System (GRS) & fusible plug type safety valve is available for Triple-Safety.
- The pressure gauge is Colour-coded Pressure Gauge to indicate the safe and danger zones.
- The repeating type safety valve is provided for better Safety and longer & error-free performance.
- The vacuum release & steam release valves are made of Stainless Steel for longer life and to avoid toxic deposits on valves due to oxidation.
- The handles is made from heat resistant, impact tolerant, finest quality bakelite, fitted with Stainless Steel screws to avoid easy loosening of handles.
- The inside tripod stand is made of Stainless Steel with extra height to avoid wetting of instruments.
- The outer stand is well designed to provide proper stability to the Autoclave Wing Nut.
- We are Manufacturer, Exporters and Suppliers of Autoclave Wing Nut Non Electric Made of Aluminium at wholesale competitive prices from India.

Model No	Size in inch	Size in mm	Capacity (Ltr)
AUWN 151	9 X 9 Inch	230 X 230 mm	10 Ltr
AUWN 152	9 X 13 Inch	230 X 330 mm	15 Ltr
AUWN 153	12 X 12 Inch	300 X 300 mm	20 Ltr
AUWN 154	12 X 12 Inch	300 X 300 mm	21 Ltr
AUWN 155	12 X 14 Inch	300 X 350 mm	25 Ltr
AUWN 156	12 X 15 Inch	300 X 375 mm	27 Ltr
AUWN 157	12 X 16 Inch	300 X 400 mm	30 Ltr
AUWN 158	14 X 15 Inch	350 X 375 mm	39 Ltr
AUWN 159	12 X 22 Inch	305 X 560 mm	40 Ltr