



Autoclave Wing Nut Stainless Steel Non Electric (MODEL NO : AUWN 251)



**Autoclave Wing Nut Non Electric Made of Stainless Steel
Manufactured by Desco Medical India**



Decide with Confidence

Autoclave Wing Nut Stainless Steel Non Electric



www.descomedicalindia.com

www.autoclavesterilizersindia.com

www.descoonline.in



Autoclave Wing Nut Stainless Steel Non Electric (MODEL NO : AUWN 251)



Autoclave Wing Nut Non Electric Made of Stainless Steel Manufactured by Desco Medical India

Specifications :-

- Autoclave Wing Nut is also known as Wing Nut Sterilizers, Portable Wing Nut Autoclaves Sterilizers
- Autoclave Wing Nut is a part of Steam Sterilizing Equipment
- Autoclave Wing Nut is Non Electric and is made of Stainless Steel
- The body and lid has Seamless Construction from Extra Thick Stainless Steel Sheet Body & Lid for more safety and life long performance.
- The Anti-Bulging & Fuel Saver Base is provided for Additional Strength, Fuel Saving and Time Saving.
- The depression at centre of the lid helps in Additional Strength.
- The Gasket is ISI Marked food grade synthetic rubber gasket for better reliability, hygiene & longer life.
- The Gasket Release System (GRS) & fusible plug type safety valve is available for Triple-Safety.
- The pressure gauge is Colour-coded Pressure Gauge to indicate the safe and danger zones.
- The repeating type safety valve is provided for better Safety and longer & error-free performance.
- The vacuum release & steam release valves are made of Stainless Steel for longer life and to avoid toxic deposits on valves due to oxidation.
- The handles is made from heat resistant, impact tolerant, finest quality bakelite, fitted with Stainless Steel screws to avoid easy loosening of handles.
- The inside tripod stand is made of Stainless Steel with extra height to avoid wetting of instruments.
- The outer stand is well designed to provide proper stability to the Autoclave Wing Nut.
- We are Manufacturer, Exporters and Suppliers of Autoclave Wing Nut Non Electric Made of Stainless Steel at wholesale competitive prices from India.

Model No	Size in inch	Size in mm	Capacity (Ltr)
AUWN 251	10 X 15 Inch	250 X 375 mm	18 Ltr
AUWN 252	12 X 12 Inch	300 X 300 mm	21 Ltr
AUWN 253	12 X 15 Inch	300 X 375 mm	27 Ltr
AUWN 254	12 X 22 Inch	300 X 550 mm	39 Ltr