



Autoclave Wing Nut Stainless Steel Semi Automatic Electric (MODEL NO : AUWN 301)



Autoclave Wing Nut Semi Automatic Electric Made of Stainless Steel Manufactured by Desco Medical India



Decide with Confidence

Autoclave Wing Nut Stainless Steel Semi Automatic Electric



www.descomedicalindia.com

www.autoclavesterilizersindia.com

www.descoonline.in



Autoclave Wing Nut Stainless Steel Semi Automatic Electric (MODEL NO : AUWN 301)



Autoclave Wing Nut Semi Automatic Electric Made of Stainless Steel Manufactured by Desco Medical India

Specifications :-

Autoclave Wing Nut Semi Automatic Electric has full Stainless Steel 304 Grade Construction.

Semi Automatic Autoclave has Seamless diepressed deep drawn chamber

µp based PID controller is fitted in Semi Automatic Autoclave with dual display for Temp. & Time

Microprocessor Digital temperature controller cum timer is fitted to the Autoclave

In case of Over heat there is an Auto shut down at 1370 C

Unique Shape heater of outer SS tubing is fitted to the Autoclave and is designed for optimum chamber usage

Food grade silicon rubber gasket for better durability

Temperature cum pressure gauge is Calibrated

Monitoring gauge is Color Coded

The Operating pressure is 1.2/2.1kg/cm²

Operating Temperature is 1210 C/1340 C

The Sterilization time is 0-99 mins

The Voltage is 230 V

Safety Valve is Stainless Steel

Vacuum Breaker are provided to remove air packets in the chamber

We are Manufacturer, Exporters and Suppliers of Autoclave Wing Nut Stainless Steel Semi Automatic Electric at wholesale competitive prices from India.

Model No	Size in inch	Size in mm	Capacity (Ltr)	Heating Element	PSI
AUWN 301	8 X 8 Inch	200 X 200 mm	6 Ltr	1 KW	15 PSI
AUWN 302	10 X 10 Inch	250 X 250 mm	12 Ltr	1.5 KW	15 PSI
AUWN 303	12 X 12 Inch	300 X 300 mm	21 Ltr	2 KW	15 PSI
AUWN 304	12 X 12 Inch	300 X 300 mm	21 Ltr S	2 KW	15 PSI
AUWN 311	8 X 8 Inch	200 X 200 mm	6 Ltr	1 KW	30 PSI
AUWN 312	10 X 10 Inch	250 X 250 mm	12 Ltr	1.5 KW	30 PSI
AUWN 313	12 X 12 Inch	300 X 300 mm	21 Ltr	2 KW	30 PSI
AUWN 314	12 X 12 Inch	300 X 300 mm	21 Ltr S	2 KW	30 PSI