

Estd. 1966



DELIVERING TO THE WORLD

Autopsy Forensic Workstation (MODEL NO : AUAW 101)



Autopsy Forensic Workstation Manufactured by Desco Medical India

Specification:-

Autopsy Workstation are made from 304 Grade Stainless Steel

Autopsy Workstation is also known as Dissection Workstation and Post Mortem Workstation

Autopsy Workstation comprises of standard table with raised edges that has faucets and drains which are used to wash away collecting fluids during the internal investigation.

These work stations are fully equipped with requisite features to provide a clean hygienic and functional platform for performing autopsy procedures.

Our Autopsy workstation are designed with dissection wing and base cabinet to provide extra storage.

Our workstations have clear hygienic bacteria free work platform space and is free from all rivets and bolts.

The dimensions of the table are 225L X 75W X 90H cm

The dimensions of Dissection wing are 155L X 75W X 97H cm

There is sufficient space for wide access below the table for sanitation and cleaning

This table comes with double bowl sink

Estd. 1966



DELIVERING TO THE WORLD

Autopsy Forensic Workstation (MODEL NO : AUAW 101)



Two piece construction for easy installation

There is wrist operated stainless steel faucet with hot and cold water flow control valves

It has reverse flow hydro aspirator with built-in vacuum breaker

This workstation is equipped with large radius inside corners for easy clean up

Autopsy Workstation has three solid stainless steel sliding body supports with wooden head rest

The Engraved measurement scale on one side for easy measurement of cadaver

We are Manufacturer, Exporters and Suppliers of Autopsy Forensic Workstation at wholesale competitive prices from India.

Model No.

AUAW 101 Autopsy Forensic Workstation