

Estd. 1966



DELIVERING TO THE WORLD

# Autopsy Saw With Vacuum Dust Collector (MODEL NO : AUDC 101)



Decide with Confidence

Estd. 1966



DELIVERING TO THE WORLD

## Autopsy Saw With Vacuum Dust Collector



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

## Autopsy Saw With Vacuum Dust Collector Manufactured by Desco Medical India

### Specification:-

This versatile and compact vacuum cleaner with eco flat folding is specially designed for the suction / vacuuming of toxic dust, finely divided silica and other hazardous airborne particles, 99.97% of all particles are filtered from the exhaust that are larger than 0.5µm. For rational working, the dust removal system comes with an outlet that is specially designed for Autopsy saw. Standby-mode takes over automatic switching on and off of the suction vacuum function. The saw can be laid on the vacuum cleaner.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)

Estd. 1966



DELIVERING TO THE WORLD

# Autopsy Saw With Vacuum Dust Collector (MODEL NO : AUDC 101)



## Technical Data of Autopsy Saw

Strong Motor with at least 13,000 RPM  
12000-21000 Oscillations/out of blade  
Motor and hand piece is integrated.  
Motor is provided with long service cord with plug.  
Hand piece with safety flange permitting firm grip and stays cool during operation.  
Easily detachable and autoclavable hand piece.  
Suitable wrench to remove blades  
CE/ISO  
220-230 V 50 Hz AC

## Accessories:

Set of 2 blades for all applications.

## Technical Data of Vacuum Dust Collector

Suction Performance : 1,300 Watts  
Air flow rate : 61 l/min  
Vacuum : 244 mbar  
Main Supply : 230 Volt, 50 Hz  
Container Volume : 35 Liter

We are Manufacturer, Exporters and Suppliers of Autopsy Saw With Vacuum Dust Collector wholesale competitive prices from India.

## Model No :

**AUDC 101**

**Autopsy Saw with Vacuum Dust Collector**