

Baby Infant Incubator With Phototherapy Manufactured by Desco Medical India



Decide with Confidence

Baby Infant Incubator with Phototherapy



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.babycareequipments.com

www.descoonline.in

Specification:-

Baby Infant Incubator has temperature controlled unit which is operated by advanced micro computer technology.

The digital display indication provides Set Temperature, Air Temperature and Skin Temperature.

Two function control modes provide flexibility while optimizing the thermal environment.

In addition the skin temperature display is larger than the set/air temperature display.

For easy cleaning & service the heater and the control unit is removable.

Acrylic hood helps full visibility & accessibility of the baby.



Baby Infant Incubator With Phototherapy (MODEL NO : BCBI 105)



A full complement of alarm provides safety convenience, with both audible and visual indicators in case of :

Probe Failure / Power Failure / High Temperature / Low Temperature

Skin Temperature display can be converted from C to °F with the help of a button.

A big LCD display is provided for more specific information & visualization of the heater power in different modes.

The High & Low Temperature alarm settings can be adjusted according to the Child condition.

Humidification is provided with the help of two water trays provided at the back of the unit.

Three side trays are provided at the bottom of the unit to keep consumables.

The unit rests on four big castors with brakes.

It comes with Phototherapy Unit.

Optional at Extra Cost :-

Oxygen Delivery System, Electric Slow Suction

We are Manufacturer, Exporters and Suppliers of Baby Infant Incubator With Phototherapy wholesale competitive prices from India

MODEL NO:

BCBI 105 BABY INFANT INCUBATOR WITH PHOTOTHERAPY CFL

BCBI 106 BABY INFANT INCUBATOR WITH PHOTOTHERAPY LED