



# Bacteriological Incubator Stainless Steel Chamber (MODEL NO : LEIN 151)



Decide with Confidence

## Bacteriological Incubator Stainless Steel Chamber

Estd. 1966



DELIVERING TO THE WORLD



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

## Bacteriological Incubator Stainless Steel Chamber Manufactured by Desco Medical

### India

#### Specification:-

Incubator has standard double wall fabrication.

Inner chamber is made out of stainless steel sheet and outer made out of thick mild steel sheet finished in white stoving enamel paint.

The gap is filled with special glass wool for temperature insulation.

Double wall door has double window for observation in the chamber.

Temperature ranges 5oC above ambient to 80oC + 2oC controlled by a thermostat.

L- shaped thermometer is provided for taking temperature readings.

Panel is provided with pilot lamps for mains & heaters with on-off switch and temperature Control knob.

Incubator works on 220 V AC 50 Hz single phase.

Heating is achieved from 3 sides.

Perforated removable shelves are provided.

Optional Accessories at Extra Cost :-

i) Air Circulating fan.

ii) Timer 0-24 hours.

iii) Digital temperature Indicator Cum Controller.

iv) Microprocessor PID digital temperature Indicator Cum Controller with timer

v) Digital Temperature Indicator

We are Manufacturer, Exporters and Suppliers of Bacteriological Incubator Stainless Steel Chamber at wholesale competitive prices from India.

Model No	Chamber Size (HXWXD) MM	Chamber Size (HXWXD) INCHES	No. of Shelves	Capacity (Approx)
LEIN 151	300 x 300 x 300 mm	12 x 12 x 12 inch	2	28 Ltrs
LEIN 152	350 x 350 x 350 mm	14 x 14 x 14 inch	2	45 Ltrs
LEIN 153	400 x 400 x 400 mm	16 x 16 x 16 inch	2	65 Ltrs
LEIN 154	450 x 450 x 450 mm	18 x 18 x 18 inch	2	95 Ltrs
LEIN 155	600 x 450 x 450 mm	24 x 18 x 18 inch	2	125 Ltrs
LEIN 156	600 x 600 x 600 mm	24 x 24 x 24 inch	2	224 Ltrs
LEIN 157	900 x 600 x 450 mm	36 x 24 x 18 inch	2	252 Ltrs
LEIN 158	900 x 600 x 600 mm	36 x 24 x 24 inch	2	336 Ltrs
LEIN 159	900 x 900 x 900 mm	36 x 36 x 36 inch	3	720 Ltrs