

Estd. 1966



DELIVERING TO THE WORLD

Blood Bank Refrigerated Centrifuge (MODEL NO : BBRC 101)



Blood Bank Refrigerated Centrifuge Manufactured by Desco Medical India



Specification:-

Features : Fully Automatic

Display : Graphic LCD Display of set and process parameters.

Motor Drive : Brush less induction motor VF drive

Memory : Programmable memory with tamper proof facility. Programmable Acceleration and Deceleration time range from 240 seconds to 480 seconds

Lid Lock : Safety lid lock with Indicator

Rotor: Swing out wind shield rotor With Metal Buckets suitable adapters for ,6 blood bags, 8 Blood Bags, 12 Blood Bags. Removable plastic cups to hold single / double blood bags with stainless steel clips

Cooling unit : Kirloskar brand ISI mark compressor, is CFC free. Condenser & Cooling units are made of copper tube & aluminium finns.

Temperature Range : 20 deg C to 30 deg C. LC: 1 deg C, Micro processor controlled rotor temperature within 1deg C.

Programmable Time : 99 Minutes with minimum resolution of 1 minute.



Blood Bank Refrigerated Centrifuge (MODEL NO : BBRC 101)



Blood Bank Refrigerated Centrifuge Manufactured by Desco Medical India

Specification:-

Technical :

Operation : Programmable type

Inner Chamber : Stainless steel easy to clean, corrosion resistant With provision for Drainage of condensed water.

Motor imbalance detection : Automatic shutdown of rotor movement with unbalance indication.

Power supply : 220/240 volts, 50 Hz Single Phase AC supply

Suitability : The equipment is suitable for operation from 0 deg to 40 deg C at 90 % RH Internal Electronic is suited for this ambient Condition.

Safety : Electronic safety device incorporated that will give audiovisual alarm in case of power failure.

RCF : Continues read out of working RPM & RCF values

Counter Balance: Counter balance & fall safe Lid with gas hinges

RPM : Max 4200

Maximum speed: Microprocessor controlled rotor speed to within 10 rpm of set value.

Acceleration/deceleration can be set manually.

We are Manufacturer, Exporters and Suppliers of Blood Bank Refrigerated Centrifuge at wholesale competitive prices from India

Model No	Capacity	No. of Bags
BBRC 101	500ml	6 blood bags
BBRC 102	500ml	8 blood bags
BBRC 103	500ml	12 blood bags