

Estd. 1966



DELIVERING TO THE WORLD

Clot Activator Tube Vacuum (MODEL NO : PPTU 501)



Clot Activator Tube Vacuum Manufactured by Desco Medical India



Specifications:-

Clot Activator Tube is Coated with micronised Silica particles which causes blood to clot rapidly. Clot Activator Tube provides enough and non polluted serum specimen for clinical test.

Cap Colour : Red

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently and centrifuge for 10 minutes after 15-30 minutes coagulation.

It is used in Blood Collection and storage for Biochemistry, Immunology and Serology tests.

We are Manufacturer, Exporters and Suppliers of Clot Activator Tube Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

Clot Activator Tube Vacuum (MODEL NO : PPTU 501)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 501	12.4 x 75 MM	2 ML	100 Pc.
PPTU 502	12.4 x 75 MM	3 ML	100 Pc.
PPTU 503	12.4 x 75 MM	4 ML	100 Pc.
PPTU 504	12.4 x 100 MM	5 ML	100 Pc.
PPTU 505	12.4 x 100 MM	6 ML	100 Pc.
PPTU 506	12.4 x 100 MM	7 ML	100 Pc.
PPTU 507	12.4 x 100 MM	8 ML	100 Pc.
PPTU 508	16 x 100 MM	9 ML	100 Pc.
PPTU 509	16 x 100 MM	10 ML	100 Pc.
PPTU 510	16 x 100 MM	11 ML	100 Pc.
PPTU 511	16 x 100 MM	12 ML	100 Pc.