

## Continuous Passive Motion Unit (Elbow) Supplied by Desco Medical India



**CE**



**D&B**  
Decide with Confidence

Estd. 1966



**DESCO**<sup>®</sup>

DELIVERING TO THE WORLD

### Continuous Passive Motion Unit (Elbow)





ELBOW PASSIVE MOTION



WRIST ROLL



SUPINATION PRONATION

[www.descomedicalindia.com](http://www.descomedicalindia.com)[www.descoonline.in](http://www.descoonline.in)

### Specification:-

Elbow CPM, a micro compute controlled unit, represents the latest breakthrough in CPM technology.

Solid state circuitry control panel with highly visual LED Digital readouts allows flexion and extension of Elbow joint through a prefixed range of motion and time.

Wrist rotation and Wrist roll movement can also be achieved.

Units is fitted on mobile Electrical jack to adjust height according to patient requirement.

Comes with necessary accessories for Supination / Pronation, Wrist Roll and Elbow passive motions.



# Continuous Passive Motion Unit (Elbow) (MODEL NO : PTCU 111)



## Continuous Passive Motion Unit (Elbow) Supplied by Desco Medical India

### Specification:-

#### Technical :

- \* Therapy Modes : Six, Mobilization Modes for continuous, Intermittent and Progressive passive movements.
- \* Rom Setting : Programmable from 00 to 120 degree for elbow movement and -90 to +90 for wrist movements.
- \* Pause Time : Programmable Extension / Flexion Paused from 00 to 99 seconds.
- \* Treatment Time : Programmable from 01 to 99 min.
- \* Inching Switch : Flexion or Extension starting position can be set according to the need of patient.
- \* Memory / Programme : 80 Pre-Set treatment patterns are stored in memory.
- \* Open Programs : 20 Programs are free to store your own patterns.
- \* Displays : Digital display of ROM setting Pause Time, Treatment Time & Graphic display of Mobilization modes.
- \* Height Adjustment : Electrically, by push button switch.
- \* Arm Setting : Right or Left Arm can be set & exercised by little adjustments
- \* Patient Safety : Patient can Stop / interrupt treatment with hand held switch.

We are Exporters and Suppliers of Continuous Passive Motion Unit (Elbow) at wholesale competitive prices from India.

#### MODEL NO:

**PTCU 111**      Continuous Passive Motion Unit (Elbow)