



Continuous Passive Motion Unit (Shoulder) (MODEL NO : PTCU 121)



Continuous Passive Motion Unit (Shoulder) Supplied by Desco Medical India



Decide with Confidence



Continuous Passive Motion Unit (Shoulder)



www.descomedicalindia.com

www.descoonline.in

Specification:-

Shoulder CPM, a micro compute controlled unit, represents the latest breakthrough in CPM technology.

Solid state circuitry control panel with highly visual LED Digital readouts allows passive motion of Shoulder joint through a prefixed range of motion and time.

Units is fitted on mobile Electrical jack to adjust height according to patient requirement.

Comes with all accessories.



Continuous Passive Motion Unit (Shoulder) (MODEL NO : PTCU 121)



Continuous Passive Motion Unit (Shoulder) Supplied by Desco Medical India

Specification:-

Technical :

- * Therapy Modes : Six, Mobilization Modes for continuous, Intermittent and Progressive passive movements.
- * Rom Setting : Programmable from 00 to 180 degree and -90 to +90 for all Shoulder movements.
- * Pause Time : Programmable Extension / Flexion Paused from 00 to 99 seconds.
- * Treatment Time : Programmable from 01 to 99 min.
- * Inching Switch : Starting position can be set according to the need of patient.
- * Memory / Programme : 80 Pre-Set treatment patterns are stored in memory.
- * Open Programs : 20 Programs are free to store your own patterns.
- * Displays : Digital display of ROM setting Pause Time, Treatment Time & Graphic display of Mobilization modes.
- * Height Adjustment : Electrically, by push button switch.
- * Arm Setting : Right or Left Arm can be set & exercised by little adjustments
- * Patient Safety : Patient can Stop / interrupt treatment with hand held switch.
- * Power Supply : 220+240V AC at 50-60Hz.

We are Exporters and Suppliers of Continuous Passive Motion Unit (Shoulder) at wholesale competitive prices from India.

MODEL NO:

PTCU 121

Continuous Passive Motion Unit (Shoulder)