

Estd. 1966



# Digital Arm Blood Pressure Monitor (MODEL NO : AXD-808)



DELIVERING TO THE WORLD



## Digital Arm Blood Pressure Monitor Supplied by Desco Medical India

### Specifications :-

- CDSCO Approved
- 2 Years Warranty
- PMMA Lens
- Blue Back Light
- USB Port with Cable
- Display : LCD Digital Display
- Display Size : L 60mm x W 50mm
- Display Backlight : Yes Blue Color
- Body Color : White and Blue
- Measurement Method : Oscillometric Method
- Measurement Range
- Pressure : 0 to 280 mmHg
- Pulse : 40 to 180/min

Estd. 1966



DELIVERING TO THE WORLD

# Digital Arm Blood Pressure Monitor (MODEL NO : AXD-808)



Accuracy / Calibration

Pressure :  $\pm 3$ mmHg or 2% of reading

Pulse : +5% of reading

Inflation : Automatic by Electric Pump

Deflation : Automatic Pressure release valve

Rapid Air Release : Automatic exhaust valve

Power Source : 4 "AA" Batteries 1.5V, USB Cable

Battery Life : Approx. 700 measurements (using new Alkaline Batteries) Automatic Power-off Function

Memory : Upto 90 measurements with date and time for each

Arm Cuff Size : 22-36 cm Medium

Functions:

WHO Indicator, Arrhythmia Checking, Average Function

Date, Time and Room Temperature Display when not in use.

Device Dimension : W 100mm x L 120mm x H 50mm

Contents with Product

1 x Bp Box

1 x Bp Monitor

1 x Bp Cuff 22-36 cm

1 x Product Manual

1 x USB Cable

1 x Storage Pouch

We are Exporters and Suppliers of Digital Arm Blood Pressure Monitor at wholesale competitive prices from India.

**Model No:**

**AXD-808** DIGITAL ARM BLOOD PRESSURE MONITOR