

Estd. 1966



DELIVERING TO THE WORLD

ECG Machine 3 Channel (MODEL NO : IEEC 105)



ECG Machine 3 Channel with Touch Screen Supplied by Desco Medical India



Specification:-

- 12 Lead ECG Synchronous Sampling System.
- Thermal Printer with 80mm print paper . synchronization print.
- 5.7" color touch screen, both push button and touch operation.
- Automatic measurement, calculation , analysis, waveform freezing.
- Support USB Flash disk & SD Memory.
- Data Transmitting to PC via RS232/USB.
- Patient's information can be input.
- Historical Data can be reviewed and printed.
- Built in Rechargeable li-ion battery ac/dc power conversion.
- Sleep mode to save energy and extend LCD Life.
- Portable Design with Hidden handle, Convenient to carry.
- Lead off & Lack off paper detection function.



ECG Machine 3 Channel (MODEL NO : IEEC 105)



7 level printing thick can be chosen.

Adapt to 110-230v, 50/60Hz AC power Supply

Features :

Touch Mode

Quick Input

Convenient Setting

Unlimited Storage

Technical Index :

CMRR : > 100 dB, with AC Filter

Input Circuit : Floating, Protection against defibrillator affect

Input Circuit Current : < 0.1 μ A

Resolution : 12bit/1000Hz

Operation Mode : Manual / Auto

Filter : AC. ECG Filter

Drift Filter : Anti-Drift System

Input Impodance : > 50M Ω

Patient Current Leakage : < 10 μ A

Specification of Fuse : T1, 6A 250V \emptyset 5-20

Recording Paper : Thermal Recording Paper

Calibrating Voltage : 1 mV \pm 3%

Voltage Tolerance : \pm 500 mV

Time Constant : >3.2s

Frequent Response : 0.05Hz – 150Hz

Noise Level : < 15 μ Vp-p

Threshold : < 20 μ V

Paper Speed : 25, 50mm/s (\pm 3%)

Sensitivity : 5, 10, 20mm/ mV

Recorder : High Resolution Thermal Printer

Battery : Li-ion Battery (11.1V , 3800mAh)

Dimension : 360mmX250mmX58mm

We are Exporters and Suppliers of ECG Machine 3 Channel with Touch Screen at wholesale competitive prices from India.

MODEL NO

IEEC 105 ECG Machine 3 Channel with Touch Screen