

Estd. 1966



DELIVERING TO THE WORLD

EDTA K2 Tube Non Vacuum (MODEL NO : PPTU 721)



EDTA K2 Tube Non Vacuum Manufactured by Desco Medical India



Decide with Confidence

EDTA K2 Tube Non Vacuum

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

EDTA K2 Tubes can effectively chelat the calcium ion of blood which can prevent the blood coagulation.

EDTA K2 Tube does not affect the count and size of the leukocyte and keep erythrocyte invariable.

Cap Colour : Lavender or Purple

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times immediately after drawing blood.

EDTA K2 Tube is used in Clinical hematology test & blood banking.

We are Manufacturer, Exporters and Suppliers of EDTA K2 Tube Non Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

EDTA K2 Tube Non Vacuum (MODEL NO : PPTU 721)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 721	12.4 x 75 MM	2 ML	100 Pc.
PPTU 722	12.4 x 75 MM	3 ML	100 Pc.
PPTU 723	12.4 x 75 MM	4 ML	100 Pc.
PPTU 724	12.4 x 100 MM	5 ML	100 Pc.
PPTU 725	12.4 x 100 MM	6 ML	100 Pc.
PPTU 726	12.4 x 100 MM	7 ML	100 Pc.
PPTU 727	12.4 x 100 MM	8 ML	100 Pc.
PPTU 728	16 x 100 MM	9 ML	100 Pc.
PPTU 729	16 x 100 MM	10 ML	100 Pc.
PPTU 730	16 x 100 MM	11 ML	100 Pc.
PPTU 731	16 x 100 MM	12 ML	100 Pc.