

Estd. 1966



DELIVERING TO THE WORLD

EDTA K2 Tube Vacuum (MODEL NO : PPTU 701)



EDTA K2 Tube Vacuum Manufactured by Desco Medical India



Decide with Confidence

EDTA K2 Tube Vacuum

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

EDTA K2 Tubes can effectively chelat the calcium ion of blood which can prevent the blood coagulation.

EDTA K2 Tube does not affect the count and size of the leukocyte and keep erythrocyte invariable.

Cap Colour : Lavender or Purple

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times immediately after drawing blood.

EDTA K2 Tube is used in Clinical hematology test & blood banking.

We are Manufacturer, Exporters and Suppliers of EDTA K2 Tube Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

EDTA K2 Tube Vacuum (MODEL NO : PPTU 701)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 701	12.4 x 75 MM	2 ML	100 Pc.
PPTU 702	12.4 x 75 MM	3 ML	100 Pc.
PPTU 703	12.4 x 75 MM	4 ML	100 Pc.
PPTU 704	12.4 x 100 MM	5 ML	100 Pc.
PPTU 705	12.4 x 100 MM	6 ML	100 Pc.
PPTU 706	12.4 x 100 MM	7 ML	100 Pc.
PPTU 707	12.4 x 100 MM	8 ML	100 Pc.
PPTU 708	16 x 100 MM	9 ML	100 Pc.
PPTU 709	16 x 100 MM	10 ML	100 Pc.
PPTU 710	16 x 100 MM	11 ML	100 Pc.
PPTU 711	16 x 100 MM	12 ML	100 Pc.