

Estd. 1966



EDTA K3 Tube Non Vacuum (MODEL NO : PPTU 761)



DELIVERING TO THE WORLD

EDTA K3 Tube Non Vacuum Manufactured by Desco Medical India



Decide with Confidence

EDTA K3 Tube Non Vacuum

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

EDTA K3 Tubes can effectively chelate the calcium ion of blood which can prevent the blood coagulation.

EDTA K3 Tube does not affect the count and size of the leukocyte and keep erythrocyte invariable.

Cap Colour : Lavender or Purple

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times immediately after drawing blood.

EDTA K3 Tube is used in Clinical hematology test & blood banking.

We are Manufacturer, Exporters and Suppliers of EDTA K3 Tube Non Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

EDTA K3 Tube Non Vacuum (MODEL NO : PPTU 761)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 761	12.4 x 75 MM	2 ML	100 Pc.
PPTU 762	12.4 x 75 MM	3 ML	100 Pc.
PPTU 763	12.4 x 75 MM	4 ML	100 Pc.
PPTU 764	12.4 x 100 MM	5 ML	100 Pc.
PPTU 765	12.4 x 100 MM	6 ML	100 Pc.
PPTU 766	12.4 x 100 MM	7 ML	100 Pc.
PPTU 767	12.4 x 100 MM	8 ML	100 Pc.
PPTU 768	16 x 100 MM	9 ML	100 Pc.
PPTU 769	16 x 100 MM	10 ML	100 Pc.
PPTU 770	16 x 100 MM	11 ML	100 Pc.
PPTU 771	16 x 100 MM	12 ML	100 Pc.