

Estd. 1966



DELIVERING TO THE WORLD

# Electric C-Arm Table with Top Slide (MODEL NO : OTCA 115)



## Electric C-Arm Table with Top Slide Manufactured by Desco Medical India



Decide with Confidence



### Electric C-Arm Table with Top Slide

Estd. 1966



DELIVERING TO THE WORLD



#### INCLUDED Standard Accessories



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

[www.operatingtablesindia.com](http://www.operatingtablesindia.com)

#### Specification:-

- \* This Electric sliding top O.T Table is suitable for all kind of surgeries, i.e. general surgery, vascular surgery, cardio, neurology, urology, gynecology, proctology, laparoscopy, traumatic surgery, plastic surgery and others.
- \* The base & column of the table is made of high quality medical stainless steel (Grade – 304).
- \* Non- reflecting surface, antibacterial & easy to clean.
- \* We provide with a top sliding of 5inch back/forward.
- \* Height adjustment, Lateral Tilt, Trendelenburg / Reverse Trendelenburg and Flex – Reflex, Chair Position and Longitudinal Sliding top are realized by Remote control.

Estd. 1966



DELIVERING TO THE WORLD

# Electric C-Arm Table with Top Slide (MODEL NO : OTCA 115)



## Technical Specs:-

Length of Table	1900 mm
Width of Table Top	520 mm
Minimum Height (without Mattress)	Min. 750 mm
Maximum Height (without Mattress)	Max. 1000 mm
Trendelenburg	25°
Reverse Trendelenburg	25°
Lateral Tilt	±20°
Back Section	±70°
Head Section	±60°
Leg Section	0°-90°
Top Sliding	5" Reverse & 5" Forward
Floor Locking System	Center Floor Locking

We are Manufacturer, Exporters and Suppliers of Electric C-Arm Table with Top Slide at wholesale competitive prices from India.

### MODEL NO:

**OTCA 115** Electric C-Arm Table with Top Slide