

Estd. 1966



DELIVERING TO THE WORLD

Electric C-Arm Table without Top Slide (MODEL NO : OTCA 110)



Electric C-Arm Table without Top Slide Manufactured by Desco Medical India



Decide with Confidence



Estd. 1966



DELIVERING TO THE WORLD

Electric C-Arm Table without Top Slide



INCLUDED Standard Accessories



www.descomedicalindia.com

www.descoonline.in

www.operatingtablesindia.com

Specification:-

- * This Electro O.T Table is suitable for all kind of surgeries, i.e. general surgery, vascular surgery, cardio, neurology, urology, gynecology, proctology, laparoscopy, traumatic surgery, plastic surgery and others.
- * The base & column of the table is made of high quality medical stainless steel (Grade – 304).
- * Non- reflecting surface, antibacterial & easy to clean. Height adjustment, Lateral Tilt, Trendelenburg / Reverse Trendelenburg and Flex – Reflex, Chair Position are realized by Remote control.

Estd. 1966



DELIVERING TO THE WORLD

Electric C-Arm Table without Top Slide (MODEL NO : OTCA 110)



Technical Specs:-

Length of Table	1900 mm
Width of Table Top	520 mm
Minimum Height (without Mattress)	Min. 750 mm
Maximum Height (without Mattress)	Max. 1000 mm
Trendelenburg	25°
Reverse Trendelenburg	25°
Lateral Tilt	±20°
Back Section	±70°
Head Section	±60°
Leg Section	0°-90°
Floor Locking System	Center Floor Locking

We are Manufacturer, Exporters and Suppliers of Electric C-Arm Table without Top Slide at wholesale competitive prices from India.

MODEL NO:

OTCA 110 Electric C-Arm Table without Top Slide