

Estd. 1966



DELIVERING TO THE WORLD

Electric OT Table C-Arm (MODEL NO : OTCA 105)



Electric OT Table C-Arm Manufactured by Desco Medical India




Decide with Confidence

Electric OT Table C-Arm

Estd. 1966



DELIVERING TO THE WORLD



INCLUDED
Standard Accessories

 Arm Support	 Shoulder Support	 Side Support
 Lateral Support	 Anesthetic Screen Rod	 Knee Crutches Goepel Type

www.descomedicalindia.com

www.descoonline.in

www.operatingtablesindia.com

Specification:-

- * We are providing SMPS based controlbox which can handle better Voltage Fluctuations, When SMPS is idle there is zero current consumption and saves power.
- * We provide remote which is ergonomically designed to have an easy and better grip with spiral cable and it has LED backlit switch with symbolic position figures making it convenient to use even in the dark.
- * Stainless Steel Fitting
- * Detachable head, leg & Pelvic section
- * Base & Cover made of attractive, easy to clean Stainless Steel.
- * Sophisticated mechanics provide smooth, step less articulation of the table Top for precise patient positioning.
- * Standby Battery Back-up with recharging circuit.
- * We provide Anti Bacterial and Anti Fungal Mattress.



Electric OT Table C-Arm (MODEL NO : OTCA 105)



Technical Specs:-

Length of Table	1900 mm
Width of Table Top	520 mm
Minimum Height (without Mattress)	Min. 750 mm
Maximum Height (without Mattress)	Max. 1000 mm
Trendelenburg	25°
Reverse Trendelenburg	25°
Lateral Tilt	±20°
Back Section	±60°
Head Section	±60°
Leg Section	0°-90°
Floor Locking System	Manual

We are Manufacturer, Exporters and Suppliers of Electric OT Table C-Arm at wholesale competitive prices from India.

MODEL NO:

OTCA 105 Electric OT Table C-Arm