

Estd. 1966



DELIVERING TO THE WORLD

Electro Hydraulic OT Table (MODEL NO : OTEH 105)



Electro Hydraulic OT Table Manufactured by Desco Medical India




Decide with Confidence

Estd. 1966



DELIVERING TO THE WORLD

Electro Hydraulic OT Table



**INCLUDED
Standard Accessories**

 Arm Support	 Shoulder Support	 Side Support
 Lateral Support	 Anesthetic Screen Rod	 Knee Crutches Goepel Type


Height Adjustment - Down


Height Adjustment - Up


Trendelenburg


Reverse Trendelenburg


Lateral Tilt - Left


Lateral Tilt - Right


Flex Position


Reflex Position

www.descomedicalindia.com

www.descoonline.in

www.operatingtablesindia.com

Specification:-

- * Stainless Steel Fitting.
- * Detachable head, leg & Pelvic section.
- * Five Sectional Radio-translucent Top.
- * Base & Cover made of attractive, easy to clean Stainless Steel.
- * Sophisticated mechanics provide smooth, step less articulation of the table Top for precise patient positioning.
- * Stand by Battery Back-up with recharging circuit.
- * This Electro sliding top O.T Table is suitable for all kind of surgeries, i.e. general surgery, vascular surgery, cardio, neurology, urology, gynecology, proctology, laparoscopy, traumatic surgery, plastic surgery and others.



DELIVERING TO THE WORLD

Electro Hydraulic OT Table (MODEL NO : OTEH 105)



- * The base & column of the table is made of high quality medical stainless steel (Grade - 304).
- * Non-reflecting surface, antibacterial & easy to clean.
- * We provide with a top sliding of 6inch back/forward.
- * Height adjustment, Lateral Tilts, Trendelenburg and Flex - Reflex, Chair Position and Longitudinal Sliding top are realized by Remote control.

Features :-

- * 350 kg Weight Capacity.
- * Center Floor Locking System.
- * 580 mm Ultra Low Height.
- * 450 mm Lifting Three Stage Column.

Technical Specs:-

Length of Table	2000 mm
Width of Table Top	520 mm
Minimum Height (without Mattress)	Min. 580 mm
Maximum Height (without Mattress)	Max. 1020 mm
Trendelenburg	25°
Reverse Trendelenburg	25°
Lateral Tilt (Left/Right)	±20°
Back Section	±85° Up / ±60° Down
Head Section	±80°
Leg Section	0°-90°
Top Sliding	300 mm

We are Manufacturer, Exporters and Suppliers of Electro Hydraulic OT Table at wholesale competitive prices from India.

MODEL NO:

OTEH 105

Electro Hydraulic OT Table