

Estd. 1966



DELIVERING TO THE WORLD

Embalming Machine Mobile (MODEL NO : AUEM 110)



Decide with Confidence

Embalming Machine Mobile

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Embalming Machine Mobile Manufactured by Desco Medical India

Specification:-

Embalming Machine is used for preserving human or animal remains by treating them with chemicals to forestall decomposition.

The intention is usually to make the deceased suitable for public or private viewing as part of the funeral ceremony or to keep them preserved for medical purposes in an anatomical laboratory.

Embalming has three goals such as sanitization, presentation, and preservation with restoration being an important additional factor in some instances.

Embalming can help preserve the body for a duration of many years.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)

Estd. 1966



DELIVERING TO THE WORLD

Embalming Machine Mobile (MODEL NO : AUEM 110)



Embalming has a very long and cross-cultural history with many cultures giving the embalming processes as a greater religious meaning.

Embalming Machine is made of stainless steel inner tank to store fluid capacity 10-15 liters

Embalming Machine is mounted on castors for easy movement

Grips is provided in Embalming Machine for lifting

I.V. stand is fixed for mounting cannula tubing and mains cable

Mains On & In use indication provided

Complete stainless steel outer body

Note : Addition Specification

- Specialized Noisless Circuit Motor to Deliver Fluid at 40-60 PSI
- Fluid Delivery rate approx 10 Ltr. per hour.
- Pulsating Delivery of Fluid
- Pressure Gauge to Monitor Pressure
- Digital Timer for Precise pre-set run

We are Manufacturer, Exporters and Suppliers of Embalming Machine Mobile at wholesale competitive prices from India

Model No.

AUEM 110 Embalming Machine Mobile