

Exercise Chair Multi Purpose Supplied by Desco Medical India



Specification:-

Multi Exercise Chair with overhead pulley system can be used for variety of exercises involving reciprocal motion for shoulder, elbow, knees & hips, Quadriceps & hamstring stretching and other resistive & assistive exercises.

The specifically designed chair which require little space and can be used in small departments where space is very limited.

Technical :

Construction : Square Tube of 30mm construction with weight pan assembly & storage tray.

Size : Floor Area, 100cm Long x 67cm Wide, the overhead bar extends 80cm beyond front portion of the frame.

Seat & Back : Foam padded, Back's angle can be adjusted for comfort & fitted with stabilising straps.

Exercise Chair Multi Purpose Supplied by Desco Medical India

Specification:-

Overhead Pulley Blocks : Two Sets, pulley blocks can be adjusted on overhead bar vertically.

Pulley arms can be rotated (360 degrees) to adjusted width.

Weight Plates: 10kg (1kg plates = 8+0.5kg plates = 4)

Finish : Powder Coated Finish frame with chrome plated overhead bar & weight carrier guide bars.

Accessories :

1) Pulley with rope - 2 Sets.

2) Pulley Block - 1.

3) Wire Handles - 2.

4) Arm Slings - 2.

5) Ankle Straps - 2.

6) Gripping Bar - 1.

7) Snap Hooks - 4.

We are Exporters and Suppliers of Exercise Chair Multi Purpose at wholesale competitive prices from India.

MODEL NO:

PTME 151 Exercise Chair Multi Purpose