

Estd. 1966



DELIVERING TO THE WORLD

# Face Shield Reusable (MODEL NO : MDFS 201)



Decide with Confidence

## Face Shields Reusable

Estd. 1966



DELIVERING TO THE WORLD



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

### Face Shield Reusable Manufactured by Desco Medical India

Face Shield Reusable is part of Personal Protective Equipment (PPE) which aims to protect the entire face from hazards such as flying objects, chemical splashes in laboratories and hospitals Face Shield Reusable ensures safety during treatment, diagnosis and quarantine of patients with the novel Coronavirus disease **COVID-19**

# Face Shield Reusable (MODEL NO : MDFS 201)



## Specification:-

- \* Face Shield Material : 2000 Micron (2MM) fully Transparent acrylic sheet.
- \* Size : 9 x 13 Inch.
- \* Good Quality 1 inch elastic headband of black colour.
- \* Both side protective film lamination with removable film to protect with scratches.
- \* Odour less & materials used to bring comfort to user.
- \* Low fogging face shield.
- \* Washable and Reusable Face Shield.
- \* Enough space to wear your spectacles and mask inside face shield.
- \* Coverage : Eyes, Nose, Mouth.

## Packaging Details :-

Each Face Shield comes in an individual packed with product labelling.

We are Manufacturers, Suppliers and Exporters of Face Shields Reusable at wholesale competitive prices from India.

## Model No

**MDFS 201**

**Face Shield Reusable**