

Estd. 1966



DELIVERING TO THE WORLD

Fluoride Tube Non Vacuum (MODEL NO : PPTU 801)



Fluoride Tube Non Vacuum Manufactured by Desco Medical India



Decide with Confidence

Fluoride Tube Non Vacuum

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

Fluoride Tube has additive EDTA potassium salt

EDTA potassium salt can chelate the calcium ion of blood prevents the blood coagulation.

Sodium Fluoride control the glycolysis.

Cap Colour : Grey

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently.

Fluoride Tube is used for the analysis of blood glucose & lactate.

We are Manufacturer, Exporters and Suppliers of Fluoride Tube Non Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

Fluoride Tube Non Vacuum (MODEL NO : PPTU 801)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 801	12.4 x 75 MM	2 ML	100 Pc.
PPTU 802	12.4 x 75 MM	3 ML	100 Pc.
PPTU 803	12.4 x 75 MM	4 ML	100 Pc.
PPTU 804	12.4 x 100 MM	5 ML	100 Pc.
PPTU 805	12.4 x 100 MM	6 ML	100 Pc.
PPTU 806	12.4 x 100 MM	7 ML	100 Pc.
PPTU 807	12.4 x 100 MM	8 ML	100 Pc.
PPTU 808	16 x 100 MM	9 ML	100 Pc.
PPTU 809	16 x 100 MM	10 ML	100 Pc.
PPTU 810	16 x 100 MM	11 ML	100 Pc.
PPTU 811	16 x 100 MM	12 ML	100 Pc.