



Horizontal Autoclave Cylindrical Deluxe (MODEL NO : AUHO 301)



**Horizontal Autoclave Cylindrical Deluxe Made of Stainless Steel
Manufactured by Desco Medical India**



Decide with Confidence

Horizontal Autoclave Cylindrical Deluxe



www.descomedicalindia.com

www.autoclavesterilizersindia.com

www.descoonline.in

**Horizontal Autoclave Cylindrical Deluxe Made of Stainless Steel Manufactured
by Desco Medical India**

Specifications :-

Horizontal Autoclave Cylindrical Deluxe is made of Stainless Steel 304 Grade
Horizontal Autoclave Cylindrical Deluxe has double walled inside and outside.
All these sheets are argon arcs welded.
The insulation is 63mm with mineral wool.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)

Estd. 1966



DELIVERING TO THE WORLD

Horizontal Autoclave Cylindrical Deluxe



(MODEL NO : AUHO 301)

These sterilizers are tested hydraulically as well as electrically for safe operation.
 The hinged door is made of thick Stainless Steel Plate.
 The radial locking system comprising of thick rectangular SS bars with safety mechanism, ensure complete risk free and leak proof sterilizing cycle.
 Automatic self closing device to prevent opening of the door while sterilizer is in operation and under pressure.
 The stainless steel steam generator unit is fitted with easily replaceable heating element and is fixed at the base of SS Frame.
 Water level indicates the water level inside the boiler.
 It is fitted with Pressure Gauge, Safety valve and steam release cock.
 Control panel of this unit consists of main switch alongwith the indicating lamp, a pressure control switch.
 Pressure can be set from 10 to 20 PSI.
 Separate valve for injecting the steam into the main chamber and releasing the steam after use are also provided.
 Fully automatic by PID Controller which has the following features:
 Horizontal Autoclave can be operated on manual mode either or fully automatic mode.
 Auto purging: Set Temp.
 Purging Stop: Set Temp.
 Auto Exhaust: Cool
 Elapse Time
 Set Temperature: Could be set
 Facility for temperature control
 Auto cycle stop
 Buzzer At cycle over
 Suitable for operation at 440 Volts 3 phase 4 wire 50 cycles supply.
 Option Features Available on Extra Cost:
 a) Thermograph to record sterilizing temp. and time
 b) Door at back also
 c) Data Logger, 8 channel with basic software and RS 485 PC communication alternatively with 21 CFR software
 We are Manufacturer, Exporters and Suppliers of Horizontal Autoclave Cylindrical Deluxe at wholesale competitive prices from India.

Model No.	Size in inch	Size in mm	Capacity (Ltr)
AUHO 301	18 x 24 inch	450 x 600 mm	100 Ltr
AUHO 311	22 x 30 inch	550 x 750 mm	160 Ltr
AUHO 321	16 x 48 inch	400 x 1200 mm	187 Ltr