

Estd. 1966



DELIVERING TO THE WORLD

Hydraulic C-Arm OT Table (MODEL NO : OTCA 101)



Hydraulic C-Arm OT Table Manufactured by Desco Medical India




Decide with Confidence

Hydraulic C-Arm OT Table

Estd. 1966



DELIVERING TO THE WORLD



INCLUDED Standard Accessories



www.descomedicalindia.com

www.descoonline.in

www.operatingtablesindia.com

Specification:-

- * Stainless Steel Fitting.
- * Detachable head, leg & Pelvic section.
- * Base & Cover made of attractive, easy to clean Stainless Steel.
- * Sophisticated mechanics provide smooth, step less articulation of the table Top for precise patient positioning.
- * Eccentric Pillar / In Built Kidney Bridge
- * Five Sectional Radio – Translucent.



Hydraulic C-Arm OT Table (MODEL NO : OTCA 101)



Technical Specs:-

Length of Table	1900 mm
Width of Table Top	520 mm
Minimum Height (without Mattress)	Min. 750 mm
Maximum Height (without Mattress)	Max. 1000 mm
Trendelenburg	25°
Reverse Trendelenburg	25°
Lateral Tilt	±20°
Back Section	±60°
Head Section	±60°
Leg Section	0°-90°
Floor Locking System	Manual

We are Manufacturer, Exporters and Suppliers of Hydraulic C-Arm OT Table at wholesale competitive prices from India.

MODEL NO:

OTCA 101 Hydraulic C-Arm OT Table