

Estd. 1996



OT TABLE WITH HEAD END CONTROL MINOR (MODEL NO : OTHY 120)



DELIVERING TO THE WORLD

OT Table With Head End Control Minor Manufactured by Desco Medical India


CE


BSCIC
ISO 9001 REGISTERED


D&B
Decide with Confidence


Estd. 1966
DESCO
DELIVERING TO THE WORLD

OT TABLE WITH HEAD END CONTROL MINOR



www.descomedicalindia.com

www.operatingtablesindia.com

Specification:-

OT Table With Head End Control Minor Manufactured by Desco Medical India

Minor Head End Control Operation Table is Hydraulic Multi-Function System with Translucent top.

This multi-purpose table is designed to provide the perfect solution to X-Ray and Minor Operations.

FEATURES :-

A Single action double cylinder hydraulic pump controlled by a foot pedal.

Five section translucent top with large perenial cut out for drainage tray.

Stainless steel base with Auto Locking System.

Concealed Gear Mechanism.

Blue Rubberized Mattress thickness 50mm.

Estd. 1996



OT TABLE WITH HEAD END CONTROL MINOR (MODEL NO : OTHY 120)



DELIVERING TO THE WORLD

SPECIFICATION :-

Table Top : L1930mm X W497mm
Up-Down Movement : 762-1067mm
Lifting Position : 762 to 990mm
Height Adjustment : 228mm
Lateral Tilt : 20 Degree both sides
Kidney Section Up : 150mm
Flex/Reflex : 90 Degree / 220 Degree
Trendelenburg & Reverse Position : 30 Degree / 25 Degree
Head Section Up/Down : 30 Degree / 90 Degree
Back Section Up : 80 Degree to 90 Degree

Standard Accessories :-

Side Support with pads – 1 pair
Shoulder Support with pads -1 pair
Knee Crutches with Assembly – 1 pair
Water proof rubber mattress.
Anesthesia Screen Rod – 1 pc
Arm Rest with pads – 1 pair

We are Manufacturer, Exporters and Suppliers of OT Table With Head End Control Minor at wholesale competitive prices from India.

MODEL NO:

OTHY 120 OT Table With Head End Control Minor