

## Hydrotherapy Tank Butterfly Shape Supplied by Desco Medical India



### Specification:-

Butterfly Shaped, Hydrotherapy Bath Pool is designed for permanent installations with direct connections to Hot & Cold water supply, Overflow and Drainage. The shape & dimensions of the tank allows various arm & leg exercises to be carried out without hitting to the sides. The recessed areas at the head and sides of the bath provide good access for the therapist when monitoring & assisting the patient.

Construction: Constructed of heavy gauge (14 & 16 gauge) stainless steel sheet, electrically Argon welded tank is shaped to Butterfly modified figure.

Overall dimension: 266cm x 196cm x 86cm high.

Inside dimension: Length to inset 218cm. Width upper part 183cm. Width of middle narrow part 89cm. Width of lower part 135cm. Depth is 56cm.

Turbine Ejectors, Aerators & Carriages: TWO, with 1/2 H.P. motors, Ejectors with elevators on ball bearing carriages.

Power Supply: 220 V AC.



# Hydrotherapy Tank Butterfly Shape (MODEL NO : PTHT 101)



## Hydrotherapy Tank Butterfly Shape Supplied by Desco Medical India

### Specification:-

Accessories: 2 water inlets, 1 digital thermometer, drain and overflow pipes with all necessary fitting.

\* Adjustable Stainless Steel Head-Rest,

\* Stainless Steel Body Stretcher and Sitting seat fitted with nylon cloth.

Finish: Bright polish finish inside the tank with special paint finish from outside.

We are Exporters and Suppliers of Hydrotherapy Tank Butterfly Shape at wholesale competitive prices from India.

### MODEL NO:

**PTHT 101**      Hydrotherapy Tank Butterfly Shape