

## Hydrotherapy Tank Rectangular Shape Supplied by Desco Medical India



### Specification:-

Hydrotherapy Bath Pool is designed in compact Rectangular design for permanent installations with direct connections to Hot & Cold water supply, Overflow and Drainage. The full body tank allows various arm & leg exercises. The rectangular shape tank provides good access for the therapist when monitoring & assisting the patient.

Construction: Constructed of heavy gauge (14 & 16gauge) stainless steel sheet, electrically Argon welded tank is in Rectangular Shape.

Inside dimension: Tank measures; 200cm long x 100cm wide x 50cm deep from inside

Turbine Ejectors, Aerators & Carriages: One, with 1/2 H.P. motor, Ejectors with elevators on ball bearing carriages.



# Hydrotherapy Tank Rectangular Shape (MODEL NO : PTHT 121)



## Hydrotherapy Tank Rectangular Shape Supplied by Desco Medical India

### Specification:-

Head-Rest: Adjustable Stainless Steel Head-Support.

Stretcher: Stainless Steel Body Stretcher fitted with nylon cloth.

Stool: Stainless Steel Stainless Sitting seat fitted with nylon cloth.

Thermometer: One digital thermometer to show the temperature of water.

Power Supply: 220 V AC.

Accessories: 2 water inlets drain and overflow pipes with all necessary fitting.

Finish: Bright polish finish inside the tank with special paint finish from outside.

We are Exporters and Suppliers of Hydrotherapy Tank Rectangular Shape at wholesale competitive prices from India.

### MODEL NO:

**PTHT 121**

Hydrotherapy Tank Rectangular Shape