

## Hydrotherapy Tank T Shape Supplied by Desco Medical India



### Specification:-

Hydrotherapy Bath Pool for Aquatic Therapy, is designed for permanent installations with direct connections to Hot & Cold Water Supply, Overflow & Drainage.

Construction : Constructed of heavy stainless steel sheet (2mm & 1.5mm thick), Electrically Argon welded tank is shaped to "T" modified figure.

Overall Dimension : 200cm X 185cm X 75cm high. Inside dimension : Length 185cm . Width upper Part 170cm . Width of Lower Part 90cm . Depth is 50cm.

Turbine : Two Turbine Ejector Aerator are powered by the vertically mounted 1/2 hp electric motor & mounted on a spring loaded turbine elevator for height adjustment.

Turbine Carriages : Two Carriages fitted with ball bearings for easy movement of turbines around the tank.



# Hydrotherapy Tank T Shape (MODEL NO : PTHT 111)



## Hydrotherapy Tank T Shape Supplied by Desco Medical India

### Specification:-

The heavy gauge stainless steel tank rests on 40cm (16") high tubular stand re-enforced with cross bars.

Power Supply : 220 V AC.

Accessories : One digital thermometer, 2 Water inlets, drain and overflow pipes with fitting, Adjustable stainless steel Head rest, Stainless steel body stretcher and sitting seat fitted with nylon cloth.

Finish : Bright Polish Finish inside the tank with special paint finish from outside.

We are Exporters and Suppliers of Hydrotherapy Tank T Shape at wholesale competitive prices from India.

### MODEL NO:

**PTHT 111**      Hydrotherapy Tank T Shape