

Infant Radiant Warmer Diamond with Drawer & Phototherapy manufactured by Desco Medical India

CE  **Decide with Confidence**

**Infant Radiant Warmer Diamond
with Drawer & Phototherapy**



DESCO
Estd. 1966
DELIVERING TO THE WORLD

www.descomedicalindia.com www.babycareequipments.com

Specification:-

Infant Radiant Warmer is Advance microprocessor controlled system with self check on the electronic circuitry Skin / Manual control. Digital Display of Baby Temperature / Set Temperature .Manual / Servo two function control modes provide the clinician more flexibility while optimizing the thermal environment. The Spacious baby bed incorporated easily drop down sides for maximum access to the baby. The Control Panel is also the heart for a completely integrated alarm system with both audible and visual Indicators in the event of

- ProbeFail
- Skin High
- Skin Low
- Power Fail



Infant Radiant Warmer Diamond With Drawer & Phototherapy (MODEL NO : BCRW 136)



Manual Mode

Facility of converting the warmer from servo mode to manual mode for non-stop working. The heat in the manual mode is locked at 50% to give more protection for the baby.

Heater Assembly

The Heater assembly located above the center of the support structure consists of radiant heater and reflector with in built halogen examination lamp for assessment of Infant color. The heater can be swivel led to either direction.

Baby Bed

The attached spacious baby bed incorporates easy drop down glasses for maximum access to the patient.

Castors

The system is fitted with 4" high quality castors having brakes in front 2 castors. antistatic castors

Storage

The unit comes fitted with two drawer system which slides out for keeping consumables.

Technical Specification

Height - 1920mm , Length - 890mm , Width - 860mm

Electrical Specifications

Power Supply - 220V/50Hz

Heater Power - 650 Watts

Fuse - 5 Amp

- Undersurface LED Phototherapy

We are Manufacturer, Exporters and Suppliers of Infant Radiant Warmer Diamond with Drawer & Phototherapy at wholesale competitive prices from India

MODEL NO:

BCRW 136 INFANT RADIANT WARMER DIAMOND WITH DRAWER & PHOTOTHERAPY CFL

BCRW 137 INFANT RADIANT WARMER DIAMOND WITH DRAWER & PHOTOTHERAPY LED