

Infant Radiant Warmer Silver without Phototherapy manufactured by Desco Medical India



Decide with Confidence


DELIVERING TO THE WORLD

Infant Radiant Warmer Silver Without Phototherapy



www.descomedicalindia.comwww.babycareequipments.comwww.descoonline.in

Specification:-

Infant Radiant Warmer Silver from Desco India has advance micro processor controlled system which comes with self check on the electronic circuit for Skin/ Manual control.

The Manual / Servo two function control modes provides more flexibility while optimizing the thermal environment.

It has Digital Display of Baby Temperature / Set Temperature.

The baby Bed is Spacious and can be incorporated easily drop down sides for maximum access to the baby.

The system is fitted with 50mm high quality castors having brakes in front 2 castors.

The Control Panel is used as complete integrated alarm system with both audible and visual Indicators in the event of

- Probe Fail
- Skin High
- Skin Low
- Power Fail



Infant Radiant Warmer Silver Without Phototherapy (MODEL NO : BCRW 106)



Manual Mode is used for converting the warmer from servo mode to manual mode for non-stop working.

The heat in the manual mode is locked at 50% to give more protection for the baby.

Heater Assembly is located above the center of the support structure consisting of radiant heater and reflector with inbuilt halogen examination lamp for assessment of Infant color. The heater can be swivelled to either direction.

Electrical Specifications

* Power Supply 220V/50Hz

* Heater Power 650 Watts

* Fuse 5 Amp

Product Dimensions (Approx)

* Height- 192cm

* Length- 89cm

* Width- 86cm

2 Trays below the baby bed for keeping consumables.

We are Manufacturer, Exporters and Suppliers of Infant Radiant Warmer Silver without Phototherapy at wholesale competitive prices from India.

MODEL NO:

BCRW 106 INFANT RADIANT WARMER SILVER WITHOUT PHOTOTHERAPY