

Estd. 1966

®



DELIVERING TOO THE WORLD

Infrared Vein Finder Multi Purpose (MODEL NO : BCVF 110)



Infrared Vein Finder Multi Purpose Manufactured by Desco Medical India



SPECIFICATIONS :

Weight	: 300 gms
Size	: 210 x 60 x 60 mm
Operating Temp	: 4°C to 40°C
Storage Temp	: -10°C to 50°C
Visualization Depth	: Up to 10 mm
Operating Distance	: 10 to 20 cm., optimum: 15 cm
Operating Modes	: Different modes based on patients and physiology
Usage on Battery	: 2-3 hours of active use
Charging Time	: 5 hours to full charge
Illumination Type	: Near Infrared
Projector Type	: LED, Reflective
Battery Type	: Built in, Rechargeable, Li-Ion, 3.7V, 5000 mAh
Charger Type	: Micro USB, 5V, 2A, 110-230 V
Portable	: Yes
Indicators	: Charging: Orange, Charged: Green

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)



Infrared Vein Finder Multi Purpose (MODEL NO : BCVF 110)



ADVANTAGES

- Easy access to veins.
- Dual Visualisation – On Skin Map and LCD Display.
- Accurate puncturing in single attempt.
- Nil or minimal trauma for patients.
- Save precious time during drug administration.
- Non-invasive, light-weight and portable.

FEATURES

- Reduced trauma to patients as they do not have to go through multiple puncture attempts.
- Increased accuracy allows a better cannulation reducing the risk of double puncture.
- Reduced time for the procedure as the inputs from Veinus Pro give better information as compared to the naked eye vision.
- Increased productivity for the medical staff allowing them to attend more patients.
- Lesser chances of infection as lesser number of attempts required.
- Non-invasive, non-contact operation, meaning the device does not need to touch the subject or surface at all avoiding any infection.
- Completely safe, with no LASER radiation.

We are Manufacturer, Exporters and Suppliers of Infrared Vein Finder Multi Purpose at wholesale competitive prices from India.

MODEL NO:

BCVF 110 INFRARED VEIN FINDER MULTI PURPOSE