

Estd. 1966



DELIVERING TO THE WORLD

LED 42 Examination Light (MODEL NO : OLEX 140)



LED 42 Examination Light Manufactured by Desco Medical India



Specification:-

Fully Aerodynamic German design for laminar flow, ClearVisionTechnology for better viewing during long hours of surgery, Shadow less light, Memory Function, High Quality Sterilizable Handle suitable for autoclave, Unbreakable Fire-Resistant Aluminium Body, Fast Disinfection Process, Protect s UV and IR rays and Improves light penetration for deep cavity surgeries

Technical Data:

Intensity	90,000 Lux
Shadowless	Yes
Numbers of LEDs	42 (Warm + Cool White)
Intensity Control	10-100 Digital Capacitive Touch (Fully Remote Controlled)
Fitting Options Available	Ceiling, Wall, Mobile Mounted
Spot Diameter	100-150 mm
Depth of Light Penetration	150 mm
Colour Temperature	3800-4800K
Colour Rendering Index	93 Ra
Focusing Preset	at 1 Metre
LED Average life	> 50,000 Hrs .
Diameter of Light	500 mm
Power Supply	220V / 50Hz AC
Power Consumption	50 Watt

Estd. 1966



DELIVERING TO THE WORLD

LED 42 Examination Light (MODEL NO : OLEX 140)



We are Manufacturer, Exporters and Suppliers of LED 42 Examination Light at wholesale competitive prices from India

MODEL NO:

OLEX 140 LED 42 Examination Light

OLEX 141 LED 42 Examination Light With Battery Backup