

Estd. 1966



DELIVERING TO THE WORLD

# LED Light Source (MODEL NO : IECM 130)



## LED Light Source Supplied by Desco Medical India



### Specification:-

LED Light Source exceeds 2200 lumens bright and uniform white light.

It is best alternative to 300W Xenon light with its high energy efficiency and maintenance free working.

### Features :

- LED Cool Light
- Bright white and powerful light
- Super long life > 50,000 hrs
- Energy efficient
- Maintenance free
- Homogeneous lighting
- Compact design

### Applications :

- Operation theater
- Surgical room

### Benefits :

- High brightness is comparable to 300W traditional xenon light source, can be used for

Estd. 1966



DELIVERING TO THE WORLD

# LED Light Source (MODEL NO : IECM 130 )



laparoscopy without any problem

- Long life, it could last 50000 hours max, 120times of xenon, so need not to replace the bulb for years
- Wide range power from 90-264V/AC; 50/60 Hz
- Perfect color temperature 5000K to 6000K, color rendering is over 70
- Does not emit light in UV or IR
- Customized connector, it accepts fiber light guides with active areas
- It saves energy for more than 50%
- Environmentally friendly
- Friendly design, easy operation

## **Standard Accessories :**

Light Source - 1 No.

Fibre Optic Cable - 1 No.

Mains Cable - 1 No.

User Manual - 1 No.

Technical :

LED 80W

Light spectral Visible light only

Light Colour Temperature 6000K (Sunlight only)

Power Supply 90- 264VAC (50/60Hz)

Power Consumption 100W

Refrigeration Forced Air Cooling

Over Temperature Protection Automatic LED System power - off

Dimmer Electronic, 10% - 100%

Fiber Optic Connector Storz, Optional

Dimensions 270 (W) x 230 (D) x 210 (H)mm

Weight 3.5 Kg

We are Exporters and Suppliers of LED Light Source at wholesale competitive prices from India.

**MODEL NO**

**IECM 130**

**LED Light Source**