

Estd. 1966



DELIVERING TO THE WORLD

# Lithium Heparin Tube Non Vacuum (MODEL NO : PPTU 861)



## Lithium Heparin Tube Non Vacuum Manufactured by Desco Medical India



Decide with Confidence

Estd. 1966



DELIVERING TO THE WORLD

## Lithium Heparin Tube Non Vacuum



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

### Specifications:-

Lithium Heparin Tube are filled with anticoagulant heparin  
Anticoagulant Heparin activates antithrombins, thus blocking the coagulation cascade and produce a plasma.

Cap Colour : Green

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently.

Lithium Heparin Tube is used for Routine & Special Biochemical Test from Plasma

We are Manufacturer, Exporters and Suppliers of Lithium Heparin Tube Non Vacuum at wholesale competitive prices from India.

Estd. 1966



# Lithium Heparin Tube Non Vacuum (MODEL NO : PPTU 861)



DELIVERING TO THE WORLD

<b>MODEL NO:</b>	<b>OUTER DIA X LENGTH</b>	<b>CAPACITY ML</b>	<b>PACK SIZE</b>
PPTU 861	12.4 x 75 MM	2 ML	100 Pc.
PPTU 862	12.4 x 75 MM	3 ML	100 Pc.
PPTU 863	12.4 x 75 MM	4 ML	100 Pc.
PPTU 864	12.4 x 100 MM	5 ML	100 Pc.
PPTU 865	12.4 x 100 MM	6 ML	100 Pc.
PPTU 866	12.4 x 100 MM	7 ML	100 Pc.
PPTU 867	12.4 x 100 MM	8 ML	100 Pc.
PPTU 868	16 x 100 MM	9 ML	100 Pc.
PPTU 869	16 x 100 MM	10 ML	100 Pc.
PPTU 870	16 x 100 MM	11 ML	100 Pc.
PPTU 871	16 x 100 MM	12 ML	100 Pc.