

Estd. 1966



Lithium Heparin Tube Vacuum (MODEL NO : PPTU 841)



DELIVERING TO THE WORLD

Lithium Heparin Tube Vacuum Manufactured by Desco Medical India



Decide with Confidence

Lithium Heparin Tube Vacuum



Estd. 1966

DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

- Lithium Heparin Tube are filled with anticoagulant heparin
- Anticoagulant Heparin activates antithrombins, thus blocking the coagulation cascade and produce a plasma.
- Cap Colour : Green
- Max. Vacuum Degree Tolerance : +10 %
- Invert the Tube for 6-8 times gently.
- Lithium Heparin Tube is used for Routine & Special Biochemical Test from Plasma
- We are Manufacturer, Exporters and Suppliers of Lithium Heparin Tube Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

Lithium Heparin Tube Vacuum (MODEL NO : PPTU 841)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 841	12.4 x 75 MM	2 ML	50 Pcs
PPTU 842	12.4 x 75 MM	3 ML	50 Pcs
PPTU 843	12.4 x 75 MM	4 ML	50 Pcs
PPTU 844	12.4 x 100 MM	5 ML	50 Pcs
PPTU 845	12.4 x 100 MM	6 ML	50 Pcs
PPTU 846	12.4 x 100 MM	7 ML	50 Pcs
PPTU 847	12.4 x 100 MM	8 ML	50 Pcs
PPTU 848	16 x 100 MM	9 ML	50 Pcs
PPTU 849	16 x 100 MM	10 ML	50 Pcs
PPTU 850	16 x 100 MM	11 ML	50 Pcs
PPTU 851	16 x 100 MM	12 ML	50 Pcs