

Estd. 1966

®



DELIVERING TO THE WORLD

# Micropipette Fully Autoclavable Fixed Volume (MODEL NO : MPFI 101)



**Micropipette Fully Autoclavable Fixed Volume Manufactured by Desco Medical India**



Decide with Confidence

## Micropipette Fully Autoclavable Fixed Volume

Estd. 1966

®



DELIVERING TO THE WORLD



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

Estd. 1966



DELIVERING TO THE WORLD

# Micropipette Fully Autoclavable Fixed Volume (MODEL NO : MPFI 101)



- \* Individual finely polished stainless steel piston ensure Accuracy
- \* Light weight and ergonomic design for user comfort.
- \* Quick click volume setting, separate tip ejector and clear digital read out.
- \* Each pipette is individually tested and supplied with a calibration certificate of quality.
- \* Autoclavable tip cone.
- \* Soft touch tip ejection Super blow out for micro volume pipetting (0.5  $\mu$ l to 20  $\mu$ l models) easy in recalibration and maintenance.

The slender, ribbed hand grip fits comfortably in either right or left hand.

We are Manufacturer Exporters and Suppliers of Micropipette Fully Autoclavable Fixed Volume at wholesale competitive prices from India

<b>MODEL NO:</b>	<b>VOLUME RANGE</b>	<b>ACCURACY %</b>	<b>PRECISION(CV%)</b>	<b>PACK SIZE</b>
MPFI 101	1 ul	2.5	1.8	1 Pc.
MPFI 102	2.5 ul	2.5	1.2	1 Pc.
MPFI 103	5 ul	1.5	0.8	1 Pc.
MPFI 104	10 ul	1.0	0.5	1 Pc.
MPFI 105	20 ul	0.5	0.3	1 Pc.
MPFI 106	25 ul	0.5	0.3	1 Pc.
MPFI 107	50 ul	0.3	0.3	1 Pc.
MPFI 108	100 ul	0.3	0.2	1 Pc.
MPFI 109	200 ul	0.3	0.2	1 Pc.
MPFI 110	250 ul	0.3	0.2	1 Pc.
MPFI 111	500 ul	0.3	0.2	1 Pc.
MPFI 112	1000 ul	0.3	0.2	1 Pc.
MPFI 113	2000 ul	0.3	0.2	1 Pc.
MPFI 114	5000 ul	0.3	0.2	1 Pc.
MPFI 115	10000 ul	0.3	0.2	1 Pc.