

Estd. 1966



DELIVERING TO THE WORLD

Micropipette Mini Fixed Volume (MODEL NO : MPMI 101)



Micropipette Mini Fixed Volume Manufactured by Desco Medical India



Decide with Confidence



Estd. 1966

DELIVERING TO THE WORLD

Micropipette Mini Fixed Volume



www.descomedicalindia.com

www.descoonline.in

Specifications:-

Micropipette Mini Fixed Volume is the most accurate, reliable, comfortable for extensive pipetting in today's research.

Research demands only highest quality when it comes to pipetting and precision.

Handypette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

* Individual finely polished stainless steel piston ensure Accuracy

* Light weight and ergonomic design for user comfort.

* Each pipette is individually tested and supplied with a calibration certificate of quality.

* Autoclavable tip cone.

* Mini micropipette is fully autoclavable at 121C & 15 psi. (Strict Quality Control)

We are Manufacturer Exporters and Suppliers of Micropipette Mini Fixed Volume at wholesale competitive prices from India

Estd. 1966

®



Micropipette Mini Fixed Volume (MODEL NO : MPMI 101)



DELIVERING TO THE WORLD

MODEL NO:	VOLUME RANGE	ACCURACY %	PRECISION(CV%)	PACK SIZE
MPMI 101	3 ul	1.5	1.0	1 Pc.
MPMI 102	5 ul	1.5	1.0	1 Pc.
MPMI 103	10 ul	1.0	1.0	1 Pc.
MPMI 104	20 ul	0.5	0.5	1 Pc.
MPMI 105	25 ul	0.5	0.5	1 Pc.
MPMI 106	50 ul	0.4	0.3	1 Pc.
MPMI 107	100 ul	0.3	0.3	1 Pc.
MPMI 108	200 ul	0.5	0.5	1 Pc.
MPMI 109	250 ul	0.5	0.5	1 Pc.
MPMI 110	500 ul	0.5	0.5	1 Pc.
MPMI 111	1000 ul	0.5	0.5	1 Pc.