

Estd. 1966

®



# Micropipette Non Autoclavable Fixed Volume (MODEL NO : MPFI 301)



DELIVERING TO THE WORLD

## Micropipette Non Autoclavable Fixed Volume Manufactured by Desco Medical India



Decide with Confidence

### Micropipette Non Autoclavable Fixed Volume

Estd. 1966

®



DELIVERING TO THE WORLD



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

#### Specifications:-

Micropipette Non Autoclavable Fixed Volume is the most accurate, reliable, comfortable for extensive pipetting in today's research.

Research demands only highest quality when it comes to pipetting and precision. Handypette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)



# Micropipette Non Autoclavable Fixed Volume (MODEL NO : MPFI 301)



- \* Individual finely polished stainless steel piston ensure Accuracy
- \* Light weight and ergonomic design for user comfort.
- \* Quick click volume setting, separate tip ejector and clear digital read out.
- \* Each pipette is individually tested and supplied with a calibration certificate of quality.
- \* Autoclavable tip cone.
- \* Soft touch tip ejection Super blow out for micro volume pipetting (0.5  $\mu$ l to 20  $\mu$ l models) easy in recalibration and maintenance.

The slender, ribbed hand grip fits comfortably in either right or left hand.

We are Manufacturer Exporters and Suppliers of Micropipette Non Autoclavable Fixed Volume at wholesale competitive prices from India

<b>MODEL NO:</b>	<b>VOLUME RANGE</b>	<b>ACCURACY %</b>	<b>PRECISION(CV%)</b>	<b>PACK SIZE</b>
MPFI 301	1 ul	2.5	1.8	1 Pc.
MPFI 302	2.5 ul	2.5	1.2	1 Pc.
MPFI 303	5 ul	1.5	0.8	1 Pc.
MPFI 304	10 ul	1.0	0.5	1 Pc.
MPFI 305	20 ul	0.5	0.3	1 Pc.
MPFI 306	25 ul	0.5	0.3	1 Pc.
MPFI 307	50 ul	0.3	0.3	1 Pc.
MPFI 308	100 ul	0.3	0.2	1 Pc.
MPFI 309	200 ul	0.3	0.2	1 Pc.
MPFI 310	250 ul	0.3	0.2	1 Pc.
MPFI 311	500 ul	0.3	0.2	1 Pc.
MPFI 312	1000 ul	0.3	0.2	1 Pc.
MPFI 313	2000 ul	0.3	0.2	1 Pc.
MPFI 314	5000 ul	0.3	0.2	1 Pc.
MPFI 315	10000 ul	0.3	0.2	1 Pc.