

Estd. 1966

®



DELIVERING TO THE WORLD

# Micropipette Semi Autoclavable Fixed Volume (MODEL NO : MPFI 201)



## Micropipette Semi Autoclavable Fixed Volume Manufactured by Desco Medical India



Decide with Confidence

### Micropipette Semi Autoclavable Fixed Volume

Estd. 1966

®



DELIVERING TO THE WORLD



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

#### Specifications:-

Micropipette Semi Autoclavable Fixed Volume is the most accurate, reliable, comfortable for extensive pipetting in today's research.

Research demands only highest quality when it comes to pipetting and precision. Handypette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)



# Micropipette Semi Autoclavable Fixed Volume (MODEL NO : MPFI 201)



- \* Individual finely polished stainless steel piston ensure Accuracy
- \* Light weight and ergonomic design for user comfort.
- \* Quick click volume setting, separate tip ejector and clear digital read out.
- \* Each pipette is individually tested and supplied with a calibration certificate of quality.
- \* Autoclavable tip cone.
- \* Soft touch tip ejection Super blow out for micro volume pipetting (0.5  $\mu$ l to 20  $\mu$ l models) easy in recalibration and maintenance.

The slender, ribbed hand grip fits comfortably in either right or left hand.

We are Manufacturer Exporters and Suppliers of Micropipette Semi Autoclavable Fixed Volume at wholesale competitive prices from India

<b>MODEL NO:</b>	<b>VOLUME RANGE</b>	<b>ACCURACY %</b>	<b>PRECISION(CV%)</b>	<b>PACK SIZE</b>
MPFI 201	1 ul	2.5	1.8	1 Pc.
MPFI 202	2.5 ul	2.5	1.2	1 Pc.
MPFI 203	5 ul	1.5	0.8	1 Pc.
MPFI 204	10 ul	1.0	0.5	1 Pc.
MPFI 205	20 ul	0.5	0.3	1 Pc.
MPFI 206	25 ul	0.5	0.3	1 Pc.
MPFI 207	50 ul	0.3	0.3	1 Pc.
MPFI 208	100 ul	0.3	0.2	1 Pc.
MPFI 209	200 ul	0.3	0.2	1 Pc.
MPFI 210	250 ul	0.3	0.2	1 Pc.
MPFI 211	500 ul	0.3	0.2	1 Pc.
MPFI 212	1000 ul	0.3	0.2	1 Pc.
MPFI 213	2000 ul	0.3	0.2	1 Pc.
MPFI 214	5000 ul	0.3	0.2	1 Pc.
MPFI 215	10000 ul	0.3	0.2	1 Pc.