

Estd. 1966

®



DELIVERING TO THE WORLD

Micropipette Semi Autoclavable Variable Volume (MODEL NO : MPVA 201)



Micropipette Semi Autoclavable Variable Volume Manufactured by Desco Medical India



Specifications:-

Micropipette Semi Autoclavable Variable Volume is the most accurate, reliable, comfortable for extensive pipetting in today's research.

Research demands only highest quality when it comes to pipetting and precision. Handypette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

- * Individual finely polished stainless steel piston ensure Accuracy
- * Light weight and ergonomic design for user comfort.
- * Quick click volume setting, separate tip ejector and clear digital read out.
- * Each pipette is individually tested and supplied with a calibration certificate of quality.
- * Autoclavable tip cone.
- * Soft touch tip ejection Super blow out for micro volume pipetting (0.5 μ l to 20 μ l models) easy in recalibration and maintenance.

Estd. 1966



DELIVERING TO THE WORLD

Micropipette Semi Autoclavable Variable Volume (MODEL NO : MPVA 201)



The slender, ribbed hand grip fits comfortably in either right or left hand.

We are Manufacturer Exporters and Suppliers of Micropipette Semi Autoclavable Variable Volume at wholesale competitive prices from India

MODEL NO:	VOLUME RANGE	ACCURACY %	PRECISION(CV%)	PACK SIZE
MPVA 201	0.2-2.0 ul	12.0	6.0	1 Pc.
		2.5	1.5	
		1.5	0.7	
MPVA 202	0.5 - 10 ul	5.0	2.8	1 Pc.
		1.5	0.8	
		1.0	0.4	
MPVA 203	2.0 - 20 ul	5.0	1.5	1 Pc.
		1.0	0.6	
		1.0	0.3	
MPVA 204	5.0 - 50 ul	2.5	6.0	1 Pc.
		0.7	1.5	
		0.7	0.7	
MPVA 205	10 - 100 ul	2.5	2.8	1 Pc.
		0.8	0.8	
		0.8	0.4	
MPVA 206	20 - 100 ul	2.5	1.5	1 Pc.
		1.0	0.6	
		0.6	0.3	
MPVA 207	100 - 1000 ul	1.6	0.3	1 Pc.
		0.6	0.2	
		0.5	0.2	
MPVA 208	200 - 2000 ul	1.5	0.3	1 Pc.
		0.6	0.2	
		0.6	0.2	
MPVA 209	0.5 - 5 ml	1.5	0.3	1 Pc.
		0.6	0.2	
		0.6	0.2	
MPVA 210	1 - 10 ml	1.5	0.3	1 Pc.
		0.6	0.2	
		0.6	0.2	